

# Time Series Analysis in



## Outline: 6+ Hours of Edification

- Philosophy (e.g., theory without equations)
- Sample FMRI data
- Theory underlying FMRI analyses: the **HRF**
- “Simple” or “Fixed Shape” regression analysis
  - ★ Theory and Hands-on examples
- “Deconvolution” or “Variable Shape” analysis
  - ★ Theory and Hands-on examples
- Advanced Topics (followed by brain meltdown)

**Goals: Conceptual Understanding + Prepare to Try It Yourself**

## Data Analysis Philosophy

- **Signal** = Measurable response to stimulus
- **Noise** = Components of measurement that interfere with detection of signal
- Statistical detection theory:
  - ★ **Understand** relationship between stimulus & signal
  - ★ Characterize noise statistically
  - ★ Can then devise methods to distinguish noise-only measurements from signal+noise measurements, and assess the methods' reliability
  - ★ Methods and usefulness depend strongly on the assumptions
    - Some methods are more “robust” against erroneous assumptions than others, but may be less sensitive

## FMRI Philosophy: Signals and Noise

- FMRI [Stimulus](#)→[Signal](#) connection and [noise statistics](#) are both complex and poorly characterized
- Result: there is no “**best**” way to analyze FMRI time series data: there are only “**reasonable**” analysis methods
- To deal with data, must make some assumptions about the signal and noise
- Assumptions will be wrong, but must do *something*
- Different kinds of experiments require different kinds of analyses
  - ★ Since signal models and questions you ask about the signal will vary
  - ★ It is important to [understand](#) what is going on, so you can select and evaluate “reasonable” analyses

## Meta-method for creating analysis methods

- Write down a mathematical model connecting stimulus (or “activation”) to signal
- Write down a statistical model for the noise
- Combine them to produce an equation for measurements given signal+noise
  - ★ Equation will have unknown parameters, which are to be estimated from the data
  - ★ N.B.: signal may have zero strength (no “activation”)
- Use statistical detection theory to produce an algorithm for processing the measurements to assess signal presence and characteristics
  - ★ e.g., least squares fit of model parameters to data

## Time Series Analysis on Voxel Data


- Most common forms of fMRI analysis involve fitting an activation+BOLD model to each voxel's time series *separately* (AKA “univariate” analysis)
  - ★ Some pre-processing steps do include inter-voxel computations; e.g.,
    - spatial smoothing to reduce noise
    - spatial registration to correct for subject motion
- Result of model fits is a set of parameters at each voxel, estimated from that voxel's data
  - ★ e.g., activation amplitude ( $\beta$ ), delay, shape
  - ★ “SPM” = statistical parametric map; e.g.,  $t$  or  $F$
- Further analysis steps operate on individual SPMs
  - ★ e.g., combining/contrasting data among subjects

## Some Features of fMRI Voxel Time Series

- fMRI only measures *changes* due to neural “activity”
  - ★ Baseline level of signal in a voxel means little or nothing about neural activity
  - ★ Also, baseline level tends to drift around slowly (100 s time scale or so; mostly from small subject motions)
- Therefore, an fMRI experiment must have at least 2 different neural conditions (“tasks” and/or “stimuli”)
  - ★ Then statistically test for differences in the MRI signal level between conditions
  - ★ Many experiments: one condition is “rest”
- Baseline is modeled separately from activation signals, and *baseline model includes “rest” periods*


## Some Sample FMRI Data Time Series

- First: Block-trial FMRI data
  - ★ “Activation” occurs over a sustained period of time (say, 10 s or longer), usually from more than one stimulation event, in rapid succession
  - ★ BOLD (hemodynamic) response accumulates from multiple close activations and is large
  - ★ BOLD response is often visible in time series
  - ★ Noise magnitude about same as BOLD response
- Next 2 slides: same brain voxel in 3 (of 9) EPI runs
  - ★ **black curve** (noisy) = data
  - ★ **red curve** (above data) = ideal model response
  - ★ **blue curve** (within data) = model fitted to data
  - ★ somatosensory task (finger being rubbed)



Block-trials: 27 s “on” / 27 s “off”; TR=2.5 s; 130 time points/run

### Same Voxel: Run 3 and Average of all 9




⇒ Activation amplitude and shape are variable! Why???


## More Sample FMRI Data Time Series

- Second: Event-related FMRI
  - ★ “Activation” occurs in single relatively brief intervals
  - ★ “Events” can be randomly or regularly spaced in time
    - If events are randomly spaced in time, signal model itself looks noise-like (to the pitiful human eye)
  - ★ BOLD response to stimulus tends to be weaker, since fewer nearby-in-time “activations” have overlapping signal changes (hemodynamic responses)
- Next slide: Visual stimulation experiment

“Active” voxel shown in next slide




## Two Voxel Time Series from Same Run



Lesson: ER-FMRI activation is not obvious via casual inspection

## Hemodynamic Response Function (HRF)

- **HRF** is the idealization of measurable FMRI signal change responding to a single activation cycle (up and down) from a stimulus in a voxel



Response to brief activation (< 1 s):

- delay of 1-2 s
- rise time of 4-5 s
- fall time of 4-6 s
- model equation:

$$h(t) \propto t^b e^{-t/c}$$


- $h(t)$  is signal change  $t$  seconds **after** activation

1 Brief Activation (Event)

## Linearity of HRF

- Multiple activation cycles in a voxel, closer in time than duration of HRF:

★ Assume that overlapping responses add




- Linearity is a pretty good assumption
- But not apparently perfect — about 90% correct
- Nevertheless, is widely taken to be true and is the basis for the “general linear model” (GLM) in FMRI analysis

3 Brief Activations

## Linearity and Extended Activation

- Extended activation, as in a block-trial experiment:

★ HRF accumulates over its duration ( $\approx 10$  s)




- **Black** curve = response to a single brief stimulus
- **Red** curve = activation intervals
- **Green** curve = summed up HRFs from activations
- Block-trials have larger BOLD signal changes than event-related experiments

2 Long Activations (Blocks)



## Convolution Signal Model

- FMRI signal we look for in each voxel is taken to be sum of the individual trial HRFs
  - ★ Stimulus timing is assumed known (or measured)
  - ★ Resulting time series (**blue** curves) are called the **convolution** of the HRF with the stimulus timing
- Must also allow for baseline and baseline drifting
  - ★ Convolution models only the FMRI signal **changes**



- Real data starts at and returns to a nonzero, slowly drifting baseline

## Simple Regression Models

- Assume a fixed shape  $h(t)$  for the HRF
  - ★ e.g.,  $h(t) = t^{8.6} \exp(-t/0.547)$  [MS Cohen, 1997]
  - ★ Convolved with stimulus timing (e.g., AFNI program **waver**), get ideal response function  $r(t)$
- Assume a form for the baseline
  - ★ e.g.,  $a + b \cdot t$  for a constant plus a linear trend
- In each voxel, fit data  $Z(t)$  to a curve of the form
 
$$Z(t) \approx a + b \cdot t + \beta \cdot r(t)$$


←

The signal model!

  - $a, b, \beta$  are unknown parameters to be calculated in each voxel
  - $a, b$  are “nuisance” parameters
  - $\beta$  is amplitude of  $r(t)$  in data = “how much” BOLD




## Simple Regression: Example



- Necessary baseline model complexity depends on duration of **continuous** imaging — e.g., 1 parameter per ~150 seconds

## Duration of Stimuli - Important Caveats

- Slow baseline drift (time scale 100 s and longer) makes doing fMRI with long duration stimuli difficult
  - Learning experiment, where the task is done continuously for ~15 minutes and the subject is scanned to find parts of the brain that adapt during this time interval
  - Pharmaceutical challenge, where the subject is given some psychoactive drug whose action plays out over 10+ minutes (e.g., cocaine, ethanol)
- Multiple very short duration stimuli that are also very close in time to each other are very hard to tell apart, since their HRFs will have 90-95% overlap
  - Binocular rivalry, where percept switches ~ 0.5 s



## Multiple Stimuli = Multiple Regressors

- Usually have more than one class of stimulus or activation in an experiment
  - ★ e.g., want to see size of “face activation” vis-à-vis “house activation”; or, “what” vs. “where” activity
- Need to model each separate class of stimulus with a separate response function  $r_1(t)$ ,  $r_2(t)$ ,  $r_3(t)$ , ....
  - ★ Each  $r_j(t)$  is based on the stimulus timing for activity in class number  $j$
  - ★ Calculate a  $\beta_j$  amplitude = amount of  $r_j(t)$  in voxel data time series  $Z(t)$
  - ★ Contrast  $\beta$ s to see which voxels have differential activation levels under different stimulus conditions
    - e.g., statistical test on the question  $\beta_1 - \beta_2 = 0$  ?


## Multiple Stimuli - Important Caveat

- You do **not** model the baseline condition
  - e.g., “rest”, visual fixation, high-low tone discrimination, or some other simple task
- fMRI can only measure **changes** in MR signal levels between tasks
  - So you need some simple-ish task to serve as a reference point
- The baseline model (e.g.,  **$a + b \cdot t$** ) takes care of the signal level to which the MR signal returns when the “active” tasks are turned off
  - Modeling the reference task explicitly would be redundant (or “collinear”, to anticipate a forthcoming jargon word)

## Multiple Stimuli - Experiment Design


- How many distinct stimuli do you need in each class? Our rough recommendations:
  - Short event-related designs: at least 25 events in each stimulus class (spread across multiple imaging runs) — and more is better
  - Block designs: at least 5 blocks in each stimulus class — 10 would be better
- While we’re on the subject: How many subjects?
  - Several independent studies agree that 20-25 subjects in each category are needed for highly reliable results
  - This number is more than has usually been the custom in fMRI-based studies!

## Multiple Regressors: Cartoon Animation



- **Red** curve = signal model for class #1
- **Green** curve = signal model for #2
- **Blue** curve =  $\beta_1 \cdot \#1 + \beta_2 \cdot \#2$  where  $\beta_1$  and  $\beta_2$  vary from 0.1 to 1.7 in the animation
- Goal of regression is to find  $\beta_1$  and  $\beta_2$  that make the blue curve best fit the data time series
- **Gray** curve =  $1.5 \cdot \#1 + 0.6 \cdot \#2 + \text{noise}$  = simulated data

## Multiple Regressors: Collinearity!!




- **Green** curve = signal model for #1
- **Red** curve = signal model for class #2
- **Blue** curve = signal model for #3
- **Purple** curve =  $\#1 + \#2 + \#3$  which is exactly = 1
- We cannot — *in principle or in practice* — distinguish sum of 3 signal models from constant baseline!!


**No** analysis can distinguish the cases  
 $Z(t) = 10 + 5 \cdot \#1$  and  
 $Z(t) = 0 + 15 \cdot \#1 + 10 \cdot \#2 + 10 \cdot \#3$   
 and an infinity of other possibilities

Collinear designs are **bad bad bad!**

## Multiple Regressors: Near Collinearity




## The Geometry of Collinearity - 1



- Trying to fit data as a sum of basis vectors that are nearly parallel doesn't work well: solutions can be huge
- Exactly parallel basis vectors would be impossible:
  - Determinant of matrix to invert would be zero

## The Geometry of Collinearity - 2



- Trying to fit data with too many regressors (basis vectors) doesn't work: no unique solution

## Equations: Notation

- Will generally follow notation of Doug Ward's manual for the AFNI program **3dDeconvolve**
- Time: continuous in reality, but in steps in the data
  - ★ Functions of continuous time are written like  $f(t)$
  - ★ Functions of discrete time expressed like  $f(\underbrace{n \cdot TR}_{=t_n})$  where  $n=0,1,2,\dots$  and  $TR$ =time step
  - ★ Usually use subscript notion  $f_n$  as shorthand
  - ★ Collection of numbers assembled in a column is a

**vector** and is printed in boldface:

$$\left\{ \begin{array}{l} \text{vector of} \\ \text{length } N \end{array} \right\} = \begin{bmatrix} f_0 \\ f_1 \\ f_2 \\ \vdots \\ f_{N-1} \end{bmatrix} = \mathbf{f}$$

$$\begin{bmatrix} A_{00} & A_{01} & \cdots & A_{0,N-1} \\ A_{10} & A_{11} & \cdots & A_{1,N-1} \\ \vdots & \vdots & \ddots & \vdots \\ A_{M-1,0} & A_{M-1,1} & \cdots & A_{M-1,N-1} \end{bmatrix} = \mathbf{A} = \{M \times N \text{ matrix}\}$$

## Equations: Single Response Function

- In each voxel, fit data  $Z_n$  to a curve of the form  

$$Z_n \approx a + b \cdot t_n + \beta \cdot r_n \quad \text{for } n=0, 1, \dots, N-1 \quad (N=\# \text{ time pts})$$
  - $a, b, \beta$  are unknown parameters to be calculated in each voxel
  - $a, b$  are “nuisance” baseline parameters
  - $\beta$  is amplitude of  $r(t)$  in data = “how much” BOLD
- Baseline model might be more complicated for long (> 150 s) continuous imaging runs:
  - $150 < T < 300$  s:  $a + b \cdot t + c \cdot t^2$
  - Longer:  $a + b \cdot t + c \cdot t^2 + \lceil T/150 \rceil$  low frequency components
  - Usually, also include as extra baseline components the estimated subject head movement time series, in order to remove residual contamination from such artifacts (will see example of this later)

≈1 param per 150 s

## Equations: Multiple Response Functions

- In each voxel, fit data  $Z_n$  to a curve of the form  

$$Z_n \approx [\text{baseline}]_n + \beta_1 \cdot r_n^{(1)} + \beta_2 \cdot r_n^{(2)} + \beta_3 \cdot r_n^{(3)} + \dots$$
  - $\beta_j$  is amplitude in data of  $r_n^{(j)} = r_j(t_n)$ ; i.e., “how much” of  $j^{\text{th}}$  response function in the data time series
  - In simple regression, each  $r_j(t)$  is derived directly from stimulus timing **and** user-chosen HRF model
    - In terms of stimulus times:  $r_n^{(j)} = \sum_{k=1}^{K_j} h(t_n - \tau_k^{(j)})$
    - If stimulus occurs on the imaging TR time-grid, stimulus can be represented as a 0-1 time series:  
 $[s_0^{(j)} \quad s_1^{(j)} \quad s_2^{(j)} \quad s_3^{(j)} \quad \dots]$  where  $s_k^{(j)} = 1$  if stimulus # $j$  is on at time  $t = k \cdot \text{TR}$ , and  $s_k^{(j)} = 0$  if # $j$  is off at that time:  

$$r_n^{(j)} = h_0 s_n^{(j)} + h_1 s_{n-1}^{(j)} + h_2 s_{n-2}^{(j)} + h_3 s_{n-3}^{(j)} + \dots = \sum_{q=0}^p h_q s_{n-q}^{(j)}$$



## Equations: Matrix-Vector Form

- Express **known** data vector as a sum of **known** columns with **unknown** coefficients:

$$\begin{bmatrix} z_0 \\ z_1 \\ z_2 \\ \vdots \\ z_{N-1} \end{bmatrix} \approx \begin{bmatrix} 1 \\ 1 \\ 1 \\ \vdots \\ 1 \end{bmatrix} \cdot a + \begin{bmatrix} 0 \\ 1 \\ 2 \\ \vdots \\ N-1 \end{bmatrix} \cdot b + \begin{bmatrix} r_0^{(1)} \\ r_1^{(1)} \\ r_2^{(1)} \\ \vdots \\ r_{N-1}^{(1)} \end{bmatrix} \cdot \beta_1 + \begin{bmatrix} r_0^{(2)} \\ r_1^{(2)} \\ r_2^{(2)} \\ \vdots \\ r_{N-1}^{(2)} \end{bmatrix} \cdot \beta_2 + \dots$$

or

$$\begin{bmatrix} z_0 \\ z_1 \\ z_2 \\ \vdots \\ z_{N-1} \end{bmatrix} \approx \begin{bmatrix} 1 & 0 & r_0^{(1)} & r_0^{(1)} & \dots \\ 1 & 1 & r_1^{(1)} & r_1^{(1)} & \dots \\ 1 & 2 & r_2^{(1)} & r_2^{(1)} & \dots \\ \vdots & \vdots & \vdots & \vdots & \ddots \\ 1 & N-1 & r_{N-1}^{(1)} & r_{N-1}^{(2)} & \dots \end{bmatrix} \begin{bmatrix} a \\ b \\ \beta_1 \\ \beta_2 \\ \vdots \end{bmatrix}$$

or

$$\mathbf{z} \approx \mathbf{R} \boldsymbol{\beta}$$

vector of data      matrix of columns      vector of coeff


the "design" matrix

' $\approx$ ' means "least squares"

$\mathbf{z}$  depends on the voxel;  $\mathbf{R}$  doesn't


## Visualizing the $\mathbf{R}$ Matrix

- Can graph columns, as shown below
  - But might have 20-50 columns
- Can plot columns on a grayscale, as shown at right
  - Easier to show many columns
  - In this plot, darker bars means larger numbers



## Solving $\mathbf{z} \approx \mathbf{R}\boldsymbol{\beta}$ for $\boldsymbol{\beta}$

- Number of equations = number of time points
  - ★ 100s per run, but perhaps 1000s per subject
- Number of unknowns usually in range 5–50
- Least squares solution:  $\hat{\boldsymbol{\beta}} = [\mathbf{R}^T \mathbf{R}]^{-1} \mathbf{R}^T \mathbf{z}$ 
  - ★  $\hat{\boldsymbol{\beta}}$  denotes an *estimate* of the true (unknown)  $\boldsymbol{\beta}$
  - ★ From  $\hat{\boldsymbol{\beta}}$ , calculate  $\hat{\mathbf{z}} = \mathbf{R}\hat{\boldsymbol{\beta}}$  as the *fitted model*




- $\mathbf{z} - \hat{\mathbf{z}}$  is the **residual time series** = noise (we hope)
- Collinearity: when matrix  $\mathbf{R}^T \mathbf{R}$  can't be inverted
  - ★ Near collinearity: when inverse exists but is huge

## Simple Regression: Recapitulation

- Choose HRF model  $h(t)$  [AKA *fixed-model regression*]
- Build model responses  $\mathbf{r}_n(t)$  to each stimulus class
  - ★ Using  $h(t)$  and the stimulus timing
- Choose baseline model time series
  - ★ Constant + linear + quadratic + movement?
- Assemble model and baseline time series into the columns of the  $\mathbf{R}$  matrix
- For each voxel time series  $\mathbf{z}$ , solve  $\mathbf{z} \approx \mathbf{R}\boldsymbol{\beta}$  for  $\hat{\boldsymbol{\beta}}$
- **Individual subject maps:** Test the coefficients in  $\hat{\boldsymbol{\beta}}$  that you care about for statistical significance
- **Group maps:** Transform the coefficients in  $\hat{\boldsymbol{\beta}}$  that you care about to Talairach space, and perform statistics on these  $\hat{\boldsymbol{\beta}}$  values

## Sample Data Analysis: Simple Regression

- Enough theory (for now: more to come later!)
- To look at the data: type `cd AFNI_data1/afni` ; then `afni`
- **Switch Underlay** to dataset `epi_r1`
  - ★ Then Sagittal **Image** and **Graph**
  - ★ **FIM**→**Pick Ideal** ; then click `afni/ideal_r1.1D` ; then **Set**
  - ★ Right-click in image, **Jump to (ijk)**, then `29 11 13`, then **Set**



- Data clearly has activity in sync with reference
- Data also has a big spike, which is very annoying
  - Subject head movement!

## Preparing Data for Analysis

- Six preparatory steps are common:
    - ★ Image registration (realignment): program `3dvolreg`
    - ★ Image smoothing: program `3dmerge`
    - ★ Image masking: program `3dClipLevel` or `3dAutomask`
    - ★ Conversion to percentile: programs `3dTstat` and `3dcalc`
    - ★ Censoring out time points that are bad: program `3dToutcount` or `3dTqual`
    - ★ Catenating multiple imaging runs into 1 big dataset: program `3dTcat`
- 
- Not all steps are necessary or desirable in any given case
  - In this first example, will only do registration, since the data obviously needs this correction

# Data Analysis Script

- In file **epi\_r1\_decon:**

```
waver -GAM
      -input epi_r1_stim.1D
      -TR 2.5
      > epi_r1_ideal.1D

3dvolreg -base 2
          -prefix epi_r1_reg
          -1Dfile epi_r1_mot.1D
          -verb
          epi_r1+orig
```

```
3dDeconvolve
  -input epi_r1_reg+orig
  -nfirst 2
  -num_stimts 1
  -stim_file 1 epi_r1_ideal.1D
  -stim_label 1 AllStim
  -tout
  -bucket epi_r1_func
  -fitts epi_r1_fitts
```

- \ • **waver** creates model time series from input stimulus timing in file **epi\_r1\_stim.1D**
- \ • Plot a 1D file to screen with **1dplot epi\_r1\_ideal.1D**
- \ • **3dvolreg** (3D image registration) will be covered in a later presentation

- \ • **3dDeconvolve** = regression code



- \ • Name of input dataset
- \ • Index of first sub-brick to process
- \ • Number of input model time series
- \ • Name of first input model time series file
- \ • Name for results in AFNI menus
- \ • Indicates to output  $t$ -statistic for  $\beta$  weights
- \ • Name of output “bucket” dataset (statistics)
- \ • Name of output model fit dataset

## Contents of .1D files

**epi\_r1\_stim.1D epi\_r1\_ideal.1D**

|     |          |
|-----|----------|
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 0   | 0        |
| 1   | 0        |
| 1   | 24.4876  |
| 1   | 122.869  |
| 1   | 156.166  |
| 1   | 160.258  |
| 1   | 160.547  |
| 1   | 160.547  |
| 1   | 160.547  |
| 1   | 160.547  |
| 1   | 160.547  |
| 0   | 136.059  |
| 0   | 37.6781  |
| 0   | 4.38121  |
| 0   | 0.288748 |
| 0   | 0        |
| 0   | 0        |
| ... | ...      |

- 1 line per time point
- TR=2.5 s
- 0=stim OFF
- 1=stim ON
- Note that “ideal” is delayed from stimulus
- Graphs at right created with **1dplot**



## To Run Script and View Results

- type **source epi\_r1\_decon** ; then wait for programs to run
- type **afni** to view what we've got
  - ★ **Switch Underlay** to epi\_r1\_reg (output from **3dvolreg**)
  - ★ **Switch Overlay** to epi\_r1\_func (output from **3dDeconvolve**)
  - ★ **Sagittal Image** and **Graph** viewers
  - ★ **FIM→Ignore→2** to have graph viewer not plot 1st 2 time pts
  - ★ **FIM→Pick Ideal** ; pick **epi\_r1\_ideal.1D** (output from **waver**)
- Define Overlay to set up functional coloring
  - **Olay→Allstim[0] Coef** (sets coloring to be from model fit  $\beta$ )
  - **Thr→Allstim[0] t-s** (sets threshold to be model fit  $t$ -statistic)
  - **See Overlay** (otherwise won't see the function!)
  - Play with threshold slider to get a meaningful activation map (e.g.,  $t=4$  is a decent threshold) — more on thresholds later

## Textual Output of the **epi\_r1\_decon** script

### • 3dvolreg output

```
++ Program 3dvolreg: AFNI version=AFNI_2005_12_30_0934 [32-bit]
++ Authored by: RW Cox
++ Reading input dataset ./epi_r1+orig.BRIK
++ Edging: x=3 y=3 z=1
++ Initializing alignment base
++ Starting final pass on 110 sub-bricks: 0..1..2..3.. *** ..106..107..108..109..
++ CPU time for realignment=8.82 s [=0.0802 s/sub-brick]
++ Min : roll=-0.086 pitch=-0.995 yaw=-0.325 dS=-0.310 dL=-0.010 dP=-0.680
++ Mean: roll=-0.019 pitch=-0.020 yaw=-0.182 dS=+0.106 dL=+0.085 dP=-0.314
++ Max : roll=+0.107 pitch=+0.090 yaw=+0.000 dS=+0.172 dL=+0.204 dP=+0.079
++ Wrote dataset to disk in ./epi_r1_reg+orig.BRIK
```

### • 3dDeconvolve output


```
++ Program 3dDeconvolve: AFNI version=AFNI_2005_12_30_0934 [32-bit]
++ Authored by: B. Douglas Ward, et al.
++ (108x3) Matrix condition [X]: 2.43095 } Quality Control: Good values
++ Matrix inverse average error = 1.3332e-14
++ Matrix setup time = 0.00 s
++ Calculations starting; elapsed time=0.502
++ voxel loop:0123456789.0123456789.0123456789.0123456789.0123456789. } Progress meter
++ Calculations finished; elapsed time=3.114
++ Wrote bucket dataset into ./epi_r1_func+orig.BRIK } Output indicators
++ Wrote 3D+time dataset into ./epi_r1_fitts+orig.BRIK
++ #Flops=4.18043e+08 Average Dot Product=4.56798
```

- If a program crashes, we'll need to see this textual output (at least)!

## More Viewing the Results

- Graph viewer: **Opt→Tran 1D→Dataset #N** to plot the model fit dataset output by **3dDeconvolve**
  - Will open the control panel for the **Dataset #N** plugin
  - Click first **Input** on ; then choose **Dataset epi\_r1\_fitts+orig**
  - Also choose **Color dk-blue** to get a pleasing plot
  - Then click on **Set+Close** (to close the plugin panel)
  - Should now see fitted time series in the graph viewer instead of data time series
  - Graph viewer: click **Opt→Double Plot→Overlay** on to make the fitted time series appear as an overlay curve
  - This tool lets you visualize the quality of the data fit
- Can also now overlay function on MP-RAGE anatomical by using **Switch Underlay** to **anat+orig** dataset
  - Probably won't want to graph the **anat+orig** dataset!

## Stimulus Correlated Movement?



- **3dvolreg** saved the motion parameters estimates into file **epi\_r1\_mot.1D**
- For fun: **1dplot epi\_r1\_mot.1D**

- Extensive “activation” (i.e., correlation of data time series with model time series) along the top of the brain is an indicator of stimulus correlated motion artifact
- Can remain even after registration, due to errors in registration, magnetic field inhomogeneities, etc.
- Can be partially removed by using the estimated movement history (from **3dvolreg**) as additional baseline model functions

## Removing Residual Motion Artifacts


- Last part of script **epi\_r1\_decon**:

```
3dDeconvolve \
  -input epi_r1_reg+orig \
  -nfirst 2 \
  -num_stimts 7 \
  -stim_file 1 epi_r1_ideal.1D \
  -stim_label 1 AllStim \
  -stim_file 2 epi_r1_mot.1D'[0]' \
  -stim_base 2 \
  -stim_file 3 epi_r1_mot.1D'[1]' \
  -stim_base 3 \
  -stim_file 4 epi_r1_mot.1D'[2]' \
  -stim_base 4 \
  -stim_file 5 epi_r1_mot.1D'[3]' \
  -stim_base 5 \
  -stim_file 6 epi_r1_mot.1D'[4]' \
  -stim_base 6 \
  -stim_file 7 epi_r1_mot.1D'[5]' \
  -stim_base 7 \
  -tout \
  -bucket epi_r1_func_mot \
  -fitts epi_r1_fitts_mot
```


These new lines add 6 regressors to the model and assign them to the baseline (-stim\_base option)

Output files: take a moment to look at results

## Some Results: Before and After



**Before:** movement parameters *are not* in baseline model



**After:** movement parameters *are* in baseline model

*t*-statistic threshold set to (uncorrected) *p*-value of  $10^{-4}$  in both images



## Setting the Threshold: Principles

- Bad things:
  - False positives — activations reported that aren't really there  $\equiv$  Type I error
  - False negatives — non-activations reported where there should be true activations found  $\equiv$  Type II error
- Usual approach in statistical testing is to control the probability of a type I error
- In FMRI, we are making many statistical tests: one per voxel ( $\approx 20,000+$ ) — the result of which is an “activation map”:
  - Voxels are colorized if they survive the thresholding process

Important Aside

## Setting the Threshold: Bonferroni

- If we set the threshold so there is a 1% chance that any given voxel is declared “active” even if its data is pure noise (FMRI jargon: “uncorrected” p-value is 0.01):
  - And assume each voxel's noise is independent of its neighbors (not really true)
  - With 30,000 voxels to threshold, would expect to get 300 false positives — this may be as many as the true activations! Situation: **Not so good.**
- Bonferroni solution: set threshold (e.g., on  $t$ -statistic) so high that uncorrected p-value is  $0.05/20000=2.5e-6$ 
  - Then have only a 5% chance that even a single false positive voxel will be reported
  - Objection: will likely lose weak areas of activation

Important Aside

## Setting the Threshold: Spatial Clustering

- Cluster-based detection lets us lower the statistical threshold and still control the false positive rate
- **Two** thresholds:
  - **First**: a per-voxel threshold that is somewhat low (so by itself leads to a lot of false positives, scattered around)
  - **Second**: form clusters of spatially contiguous (neighboring) voxels that survive the first threshold, and keep only those clusters above a volume threshold — e.g., we don't keep isolated "active" voxels
- Usually: choose volume threshold, then calculate voxel-wise statistic threshold to get the overall "corrected" p-value you want (typically, corrected  $p=0.05$ )
  - No easy formulas for this, so must use simulation: AFNI program **AlphaSim**

**Important Aside**

## AlphaSim: Clustering Thresholds

- Simulated for brain mask of 18,465 voxels
- Look for smallest cluster with corrected  $p < 0.05$

| Uncorrected<br>p-value<br>(per voxel) | Cluster Size<br>/ Corrected p<br>(uncorrelated) | Cluster Size<br>/ Corrected p<br>(correlated 5 mm) |
|---------------------------------------|---|--|
| 0.0002                                | 2 / 0.001                                       | 3 / 0.004  |
| 0.0004                                | 2 / 0.008                                       | 4 / 0.012  |
| 0.0007                                | 2 / 0.026                                       | 3 / 0.031  |
| 0.0010                                | 3 / 0.001                                       | 4 / 0.007  |
| 0.0020                                | 3 / 0.003                                       | 4 / 0.032  |
| 0.0030                                | 3 / 0.008                                       | 5 / 0.013  |
| 0.0040                                | 3 / 0.018                                       | 5 / 0.029  |
| 0.0050                                | 3 / 0.030                                       | 6 / 0.012  |
| 0.0060                                | 4 / 0.003                                       | 6 / 0.023  |
| 0.0070                                | 4 / 0.004                                       | 6 / 0.036  |
| 0.0080                                | 4 / 0.006                                       | 7 / 0.016  |
| 0.0090                                | 4 / 0.010                                       | 7 / 0.027  |
| 0.0100                                | 4 / 0.015                                       | 7 / 0.042  |

**Corresponds  
to sample data**

Can make  
activation maps  
for display with  
cluster editing  
using **3dmerge**  
program or in  
AFNI GUI (new:  
Sep 2006)

**Important Aside**

## Multiple Stimulus Classes

- The experiment analyzed here in fact is more complicated
  - ★ There are 4 related visual stimulus types
  - ★ One goal is to find areas that are differentially activated between these different types of stimuli
  - ★ We have 4 imaging runs, 108 useful time points each (skipping first 2 in each run) that we will analyze together
    - Already registered and put together into dataset **rall\_vr+orig**
  - ★ Stimulus timing files are in subdirectory **stim\_files/**
  - ★ Script file **waver\_ht2** will create HRF models for regression:
 

```
cd stim_files
waver -dt 2.5 -GAM -input scan1to4a.1D > scan1to4a_hrf.1D
waver -dt 2.5 -GAM -input scan1to4t.1D > scan1to4t_hrf.1D
waver -dt 2.5 -GAM -input scan1to4h.1D > scan1to4h_hrf.1D
waver -dt 2.5 -GAM -input scan1to4l.1D > scan1to4l_hrf.1D
cd ..
```
  - ★ Type **source waver\_ht2** to run this script
    - Might also use **1dplot** to check if input .1D files are reasonable

## Regression with Multiple Model Files

- Script file **decon\_ht2** does the job:

```
3dDeconvolve -xout -input rall_vr+orig \
-num_stimts 4 \
-stim_file 1 stim_files/scan1to4a_hrf.1D -stim_label 1 Actions \
-stim_file 2 stim_files/scan1to4t_hrf.1D -stim_label 2 Tool \
-stim_file 3 stim_files/scan1to4h_hrf.1D -stim_label 3 HighC \
-stim_file 4 stim_files/scan1to4l_hrf.1D -stim_label 4 LowC \
-concat contrasts/runs.1D \
-glt 1 contrasts/contr_AvsT.txt -glt_label 1 AvsT \
-glt 1 contrasts/contr_HvsL.txt -glt_label 2 HvsL \
-glt 1 contrasts/contr_ATvsHL.txt -glt_label 3 ATvsHL \
-full_first -fout -tout \
-bucket func_ht2
```


- Run this script by typing **source decon\_ht2** (takes a few minutes)

- Stim #1 = visual presentation of active movements
- Stim #2 = visual presentation of simple (tool-like) movements
- Stims #3 and #4 = high and low contrast gratings

## Regressors for This Script



via **1dgrayplot**  
or **-xjpeg** option



via **1dplot**  
on **-stim\_file** inputs

## Extra Features of 3dDeconvolve - 1

**-concat contrasts/runs.1D** = file that indicates where new imaging runs start

0  
108  
216  
324

**-full\_first** = put **full model** statistic first in output file, not last

**-fout -tout** = output both *F*- and *t*-statistics

- The full model statistic is an *F*-statistic that shows how well the sum of all 4 input model time series fits voxel time series data
- The individual models also will get individual *F*- and *t*-statistics indicating the significance of their individual contributions to the time series fit
  - ★ i.e.,  $F_{\text{Actions}}$  tells if model (**Actions+HighC+LowC+Tools+baseline**) explains more of the data variability than model (**HighC+LowC+Tools+baseline**) — with **Actions** omitted

## Extra Features of 3dDeconvolve - 2


```
-glt 1 contrasts/contr_AvsT.txt    -glt_label 1 AvsT
-glt 1 contrasts/contr_HvsL.txt    -glt_label 2 HvsL
-glt 1 contrasts/contr_ATvsHL.txt  -glt_label 3 ATvsHL
```

- **GLTs** are General Linear Tests
- **3dDeconvolve** provides tests for each regressor separately, but if you want to test combinations or contrasts of the  $\beta$  weights in each voxel, you need the **-glt** option
- File **contrasts/contr\_AvsT.txt** = **000000001-100**  
(one line with 12 numbers)  
8 zeros: could also write 8@0
- Goal is to test a linear combination of the  $\beta$  weights
  - ★ In this data, we have 12  $\beta$  weights: 8 baseline parameters (2 per imaging run), which are first in the  $\beta$  vector, and 4 regressor magnitudes, which are from **-stim\_file** options
  - ★ This particular test contrasts the Actions and Tool  $\beta$ s
    - tests if  $\beta_{\text{Actions}} - \beta_{\text{Tool}} \neq 0$

## Extra Features of 3dDeconvolve - 3

- File **contrasts/contr\_HvsL.txt** = **00000000001-1**
  - Goal is to test if  $\beta_{\text{HighC}} - \beta_{\text{LowC}} \neq 0$
- File **contrasts/contr\_ATvsHL.txt** = **0000000011-1-1**
  - Goal is to test if  $(\beta_{\text{Actions}} + \beta_{\text{Tool}}) - (\beta_{\text{HighC}} + \beta_{\text{LowC}}) \neq 0$
  - Regions where this statistic is significant will have had different amounts of BOLD signal change in the activity viewing tasks versus the grating viewing tasks
    - This is a way to factor out primary visual cortex
- **-glt\_label 3 ATvsHL** option is used to attach a meaningful label to the resulting statistics sub-bricks

## Results of **decon\_ht2** Script



## Statistics from 3dDeconvolve

- An *F*-statistic measures significance of how much a model component reduced the variance of the time series data
- Full *F* measures how much the signal regressors reduced the variance over just the baseline regressors (**sub-brick #0 below**)
- Individual partial-model *F*s measures how much each individual signal regressor reduced data variance over the full model with that regressor excluded (**sub-bricks #19, #22, #25, and #28 below**)
- The **Coef** sub-bricks are the  $\beta$  weights (e.g., **#17, #20, #23, #26**)
- A *t*-statistic sub-brick measure impact of one coefficient


|                    |                    |                    |
|--------------------|--------------------|--------------------|
| # 0 Full F-stat    | #13 Run #4 t^0 Coe | #26 LowC[0] Coef   |
| # 1 Run #1 t^0 Coe | #14 Run #4 t^0 t-s | #27 LowC[0] t-st   |
| # 2 Run #1 t^0 t-s | #15 Run #4 t^1 Coe | #28 LowC F-stat    |
| # 3 Run #1 t^1 Coe | #16 Run #4 t^1 t-s | #29 AvsT LC[0] coe |
| # 4 Run #1 t^1 t-s | #17 Actions[0] Coe | #30 AvsT LC[0] t-s |
| # 5 Run #2 t^0 Coe | #18 Actions[0] t-s | #31 AvsT F-stat    |
| # 6 Run #2 t^0 t-s | #19 Actions F-stat | #32 HvsL LC[0] coe |
| # 7 Run #2 t^1 Coe | #20 Tool[0] Coef   | #33 HvsL LC[0] t-s |
| # 8 Run #2 t^1 t-s | #21 Tool[0] t-st   | #34 HvsL F-stat    |
| # 9 Run #3 t^0 Coe | #22 Tool F-stat    | #35 ATvsHL LC[0] c |
| #10 Run #3 t^0 t-s | #23 HighC[0] Coef  | #36 ATvsHL LC[0] t |
| #11 Run #3 t^1 Coe | #24 HighC[0] t-st  | #37 ATvsHL F-stat  |
| #12 Run #3 t^1 t-s | #25 HighC F-stat   |                    |

## Alternative Way to Run **waver**

Instead of giving stimulus timing on TR-grid as set of 0s and 1s

- Can give the actual stimulus times (in seconds) using the **-tstim** option

★ **waver -dt 1.0 -GAM -tstim 3 12 17 | ldplot -stdin**



- If times are in a file, can use **-tstim `cat filename`** to place them on the command line after **-tstim** option

★ This is most useful for event-related experiments, where the timing of stimuli is usually given explicitly

**Note** backward single quotes

## Alternative Way to Run **3dDeconvolve**

Instead of giving stimulus timing to **waver**

- Can give the actual stimulus times (in seconds) directly to **3dDeconvolve** using the **-stim\_times** option (instead of **-stim\_file** as before)
- The program will do the equivalent of **waver** inside itself to generate the necessary column(s) in the **R** matrix
- More information in the latter part of this presentation
  - ★ Is coupled with the ideas needed for “deconvolution”
  - ★ Besides input file with stimulus times, must also specify the HRF model to be used with those times
    - That is, which shape(s) are to be placed down at each stimulus time to model the ideal response



## Deconvolution Signal Models

- Simple or Fixed-shape regression (previous):
  - ★ We fixed the shape of the HRF — amplitude varies
  - ★ Used waver to generate the signal model from the stimulus timing (or could use 3dDeconvolve directly)
  - ★ Found the amplitude of the signal model in each voxel — solution to the set linear equations =  $\beta$  weights
- Deconvolution or Variable-shape regression (next):
  - ★ We allow the shape of the HRF to vary in each voxel, for each stimulus class
  - ★ Appropriate when you don't want to over-constrain the solution by assuming an HRF shape
  - ★ **Caveat:** need to have enough time points during the HRF in order to resolve its shape

## Deconvolution: Pros and Cons

- + Letting HRF shape varies allows for subject and regional variability in hemodynamics
- + Can test HRF estimate for different shapes; e.g., are later time points more “active” than earlier?
- Need to estimate more parameters for each stimulus class than a fixed-shape model (e.g., 4-15 vs. 1 parameter=amplitude of HRF)
- Which means you need more data to get the same statistical power (assuming that the fixed-shape model you would otherwise use was in fact “correct”)
- Freedom to get any shape in HRF results can give weird shapes that are difficult to interpret

## Expressing HRF via Regression Unknowns

- The tool for expressing an unknown function as a finite set of numbers that can be fit via linear regression is an **expansion in basis functions**

$$h(t) = \beta_0 \psi_0(t) + \beta_1 \psi_1(t) + \beta_2 \psi_2(t) + \dots = \sum_{q=0}^{q=p} \beta_q \psi_q(t)$$

- ★ The basis functions  $\psi_q(t)$  are known, as is the expansion order  $p$
- ★ The unknowns to be found (in each voxel) comprises the set of weights  $\beta_q$  for each  $\psi_q(t)$
- Since  $\beta$  weights appear only by multiplying known values, and HRF only appears in final signal model by linear convolution, resulting signal model is still solvable by linear regression

## Basis Function: “Sticks”

- The set of basis functions you use determines the range of possible HRFs that you can compute
- “Stick” (or Dirac delta) functions are very flexible

★ But they come with a strict limitation

- $\delta(t)$  is 1 at  $t=0$  and is 0 at all other values of  $t$

- $\psi_q(t) = \delta(t - q \cdot \text{TR})$  for  $q=0, 1, 2, \dots, p$

$$\Rightarrow h(0) = \beta_0$$


$$\Rightarrow h(\text{TR}) = \beta_1$$

$$\Rightarrow h(2 \cdot \text{TR}) = \beta_2$$

$$\Rightarrow h(3 \cdot \text{TR}) = \beta_3$$

$$\Rightarrow \text{et cetera}$$

→  $\Rightarrow h(t) = 0$  for any  $t$  not on the TR grid



## Sticks: Good Points


- Can represent arbitrary shapes of the HRF, up and down, with ease
- Meaning of each  $\beta_q$  is completely obvious
  - ★ Value of HRF at time lag  $q \cdot TR$  after activation
- **3dDeconvolve** is set up to deal with stick functions for representing HRF, so using them is very easy
  - What is called  $p$  here is given by command line option `-stim_maxlag` in the program
  - When choosing  $p$ , rule is to estimate longest duration of neural activation after stimulus onset, then add 10-12 seconds to allow for slowness of hemodynamic response

## Sticks and TR-locked Stimuli

- **$h(t) = 0$  for any  $t$  not on the TR grid**
- This limitation means that, using stick functions as our basis set, we can only model stimuli that are “locked” to the TR grid
  - ★ That is, stimuli/activations don’t occur at fully general times, but only occur at integer multiples of TR
- For example, suppose an activation is at  $t = 1.7 \cdot TR$ 
  - ★ We need to model the response at later times, such as  $2 \cdot TR$ ,  $3 \cdot TR$ , etc., so need to model  $h(t)$  at times such as  $t = (2 - 1.7) \cdot TR = 0.3 \cdot TR$ ,  $t = 1.3 \cdot TR$ , etc., after the stimulus
- But the stick function model doesn’t allow for such intermediate times
  - **or**, can allow  $\Delta t$  for sticks to be a fraction of TR for data
  - e.g.,  $\Delta t = TR/2$ , which implies twice as many  $\beta_q$  parameters to cover the same time interval (time interval needed is set by hemodynamics)
  - then would allow stimuli that occur on TR-grid or halfway in-between

## Deconvolution and Collinearity

- Regular stimulus timing can lead to collinearity!




## 3dDeconvolve with Stick Functions

- Instead of inputting a signal model time series (e.g., created with **waver** and stimulus timing), you input the stimulus timing directly
  - ★ Format: a text file with 0s and 1s, 0 at TR-grid times with no stimulus, 1 at time with stimulus
- Must specify the maximum lag (in units of TR) that we expect HRF to last after each stimulus
  - ★ This requires you to make a judgment about the activation — brief or long?
- 3dDeconvolve** returns estimated values for each  $\beta_q$ , for each stimulus class
  - ★ Usually then use a GLT to test the HRF (or pieces of it) for significance

## Extra Features of 3dDeconvolve - 4

- `-stim_maxlag k p` = option to set the maximum lag to `p` for stimulus timing file `#k` for `k=0,1,2,...`
  - ★ Stimulus timing file input using command line option `-stim_file k filename` as before
  - ★ Can also use `-stim_minlag k m` option to set the minimum lag if you want a value `m` different from `0`
  - ★ In which case there are `p-m+1` parameters in this HRF
- `-stim_nptr k r` = option to specify that there are `r` stimulus subintervals per TR, rather than just 1
  - ★ This feature can be used to get a finer grained HRF, at the cost of adding more parameters that need to be estimated
  - Need to make sure that the input stimulus timing file (from `-stim_file`) has `r` entries per TR
  - TR for `-stim_file` and for output HRF is data  $TR \div r$

## Script for Deconvolution - The Data

- `cd AFNI_data2`
  - ★ data is in `ED/` subdirectory (10 runs of 136 images each; TR=2 s)
  - ★ script in file `@s1.analyze_ht05` (in `AFNI_data2` directory)
    - stimuli timing and GLT contrast files in `misc_files/`
  - ★ start script `now` by typing `source @s1.analyze_ht05`
    - will discuss details of script while it runs (20+ min?)
- Event-related study from Mike Beauchamp 
  - ★ 10 runs with four classes of stimuli (short videos)
    - Tools moving (e.g., a hammer pounding) - **TM**
    - People moving (e.g., jumping jacks) - **HM**
    - Points outlining tools moving (no objects, just points) - **TP**
    - Points outlining people moving - **HP**
  - ★ Goal is to find if there is an area that distinguishes natural motions (HM and HP) from simpler rigid motions (TM and TP)

- Formerly LBC/NIMH
- Now UT Houston

## Script for Deconvolution - Outline

- Examine each imaging run for outliers: **3dToutcount**
- Time shift each run's slices to a common origin: **3dTshift**
- Registration of each imaging run: **3dvolreg**
- Smooth each volume in space (136 sub-bricks per run): **3dmerge**
- Create a brain mask: **3dAutomask** and **3dcalc**
- Rescale each voxel time series in each imaging run so that its average through time is 100: **3dTstat** and **3dcalc**
  - ★ If baseline is 100, then a  $\beta_q$  of 5 (say) indicates a 5% signal change in that voxel at time lag  $q$  after stimulus
- Catenate all imaging runs together into one big dataset (1360 time points): **3dTcat**
- Compute HRFs and statistics: **3dDeconvolve**
  - ★ Each HRF will have 15 time points (lags from 0 to 14) with TR=1.0 s, since input data has TR=2.0 s and we use **-stim\_nptr k r** option with **r=2**
- Average together all points of each separate HRF to get average % change in each voxel over 14 s interval: **3dTstat**

## Script for Deconvolution - 1

```
#!/bin/tcsh
if ( $#argv > 0 ) then
    set subjects = ( $argv )
else
    set subjects = ED
endif
```

This script is designed to run analyses on a lot of subjects at once. We will only analyze the ED data here. The other subjects will be included in the Group Analysis presentation.

```
#####
# Above command will run script for all our subjects - ED, EE, EF - one after
# the other if, when we execute the script, we type: ./@s1.analyze_ht05 ED EE EF.
# If we type ./@s1.analyze_ht05 or tcsh @s1.analyze_ht05, the script runs only
# for subject ED. The user will then have to go back and edit the script so
# that 'set subjects' = EE and then EF, and then run the script for each subj.
#####
```

```
foreach subj ($subjects)
```

Loop over all subjects  
(next 2 slides)

```
cd $subj
```

First step is to change to the directory that has this subject's data

## Script for Deconvolution - 2

```
#####  
# time shift, volume register and spatially blur our datasets,  
# and remove the first two time points from each run  
#####  
  
set runs = ( `count -digits 2 1 10` )  
foreach run ( $runs )
```

Loop over imaging runs 1..10  
(loop continues on next slide)

```
    set dset = ${subj}_r${run}+orig.HEAD
```

Shorthand for dataset

```
    3dToutcount -automask ${dset} \
                > toutcount_r${run}.1D
```

Outlier check:  
By itself, **3dToutcount**  
doesn't change data!  
To plot "outlier-ness":  
**1dplot toutc\_r1.1D**

```
    3dTshift -tzero 0 -heptic \
             -prefix ${subj}_r${run}_ts \
             ${dset}
```

Interpolate each voxel's  
time series to start at the  
time of slice #0

## Script for Deconvolution - 3

```
3dvolreg -verbose \
        -base ${subj}_r01_ts+orig'[2]' \
        -prefix ${subj}_r${run}_vr \
        -1Dfile dfile.r${run}.1D \
        ${subj}_r${run}_ts+orig'[2..137]'
```

Image registration  
of each run to its  
#2 sub-brick

```
3dmerge -lblur_fwhm 4 \
        -doall \
        -prefix ${subj}_r${run}_vr_bl \
        ${subj}_r${run}_vr+orig
```

Lightly blur each 3D  
volume in each dataset  
to reduce noise and  
increase functional  
overlap among runs  
and among subjects

```
3DAutomask -dilate 1 \
        -prefix mask_r${run} \
        ${subj}_r${run}_vr_bl+orig
```

Make an "inside-the-brain"  
mask for this dataset

end → End of loop over imaging runs.  
At this point, dataset **\${subj}\_r\${run}\_vr\_bl**  
contains the data for subject **\${subj}** and imaging  
run **\${run}**, which has been time-shifted, realigned,  
and blurred; also, a brain-only mask has been made



## Script for Deconvolution - 4

```
#####
# create a union mask from those of the individual runs
#####
3dcalc -a mask_r01+orig -b mask_r02+orig -c mask_r03+orig \
      -d mask_r04+orig -e mask_r05+orig -f mask_r06+orig \
      -g mask_r07+orig -h mask_r08+orig -i mask_r09+orig \
      -j mask_r10+orig \
      -expr 'or(a+b+c+d+e+f+g+h+i+j)' \
      -prefix full_mask
```

This mask dataset will be 1 inside the largest contiguous high intensity EPI region, and 0 outside that region — this makes a brain mask

**3dcalc** program = voxel-wise “calculator” for datasets.  
 Input is 10 individual run dataset masks (1 in brain, 0 outside).  
 Output is mask which is

- 1 wherever *any* individual mask is 1,
- 0 wherever *all* individual masks are 0

## Script for Deconvolution - 5

```
#####
# - re-scale each run's mean to 100
# - use full_mask to zero out non-brain voxels
#
# If the mean is 100, and the result of 3dcalc at a voxel is 106 (at
# some time point), then one can say that voxel shows a 6% increase in
# signal activity, relative to the mean.
#####

foreach run ( $runs )

    3dTstat -prefix mean_r${run} \
           ${subj}_r${run}_vr_b1+orig

    3dcalc -a ${subj}_r${run}_vr_b1+orig \
           -b mean_r${run}+orig \
           -c full_mask+orig \
           -expr "(a/b * 100) * c" \
           -prefix scaled_r${run}

    rm -f mean_r${run}+orig*

end
```

Mean of the run<sup>th</sup> dataset,  
 through time: run=1..10

- Divide each voxel value ('a') by its temporal mean ('b') and multiply by 100
- Result will have temporal mean of 100
- Voxels not in the mask will be set to 0 (by 'c')

## Script for Deconvolution - 6

```
3dTcat -prefix ${subj}_all_runs \
      scaled_r??+orig.HEAD
```

“Gluing” the runs together, since **3dDeconvolve** only operates on one input dataset at a time

```
cat dfile.r??1D > dfile.all1D
```

Also “glue” together the movement parameters output from **3dvolreg**

```
#=====
# move unloved run data into separate directories
#=====
```

```
mkdir runs_orig runs_temp
```

```
mv ${subj}_r*_vr* ${subj}_r*_ts* scaled* \
    dfile.r??1D toutcount* runs_temp
```

Gets this stuff out of the way so that we don't see it when we run AFNI later

```
mv ${subj}_r* runs_orig
```

## Script for Deconvolution - 7

```
3dDeconvolve -polort 2 \
  -input ${subj}_all_runs+orig -num_stimts 10 Input dataset \
  -concat ../misc_files/runs.1D \
  -stim_file 1 ../misc_files/all_stims.1D'[0]' 0-1 stim file #1 \
    -stim_label 1 ToolMovie \
    -stim_minlag 1 0 -stim_maxlag 1 14 -stim_nptr 1 2 \
  -stim_file 2 ../misc_files/all_stims.1D'[1]' 0-1 stim file #2 \
    -stim_label 2 HumanMovie \
    -stim_minlag 2 0 -stim_maxlag 2 14 -stim_nptr 2 2 \
  -stim_file 3 ../misc_files/all_stims.1D'[2]' 0-1 stim file #3 \
    -stim_label 3 ToolPoint \
    -stim_minlag 3 0 -stim_maxlag 3 14 -stim_nptr 3 2 \
  -stim_file 4 ../misc_files/all_stims.1D'[3]' 0-1 stim file #4 \
    -stim_label 4 HumanPoint \
    -stim_minlag 4 0 -stim_maxlag 4 14 -stim_nptr 4 2 \
```

- 4 time series models: one for each the 4 different classes of events
- All stimuli time series in 1 file with 4 columns: **../misc\_files/all\_stims.1D**
  - Selectors like **'[2]'** pick out a particular column
  - Each stimulus input and HRF output is sampled at  $TR/2 = 1.0$  s
    - Due to the use of **-stim\_nptr k 2** for each **k**
  - Lag from 0 to 14 s is about right for HRF to a brief stimulus
- **-stim\_label** option: names used in AFNI and below in **-gltsym** options

## Script for Deconvolution - 8

```

-stim_file 5  dfile.all.1D'[0]' -stim_base 5
-stim_file 6  dfile.all.1D'[1]' -stim_base 6
-stim_file 7  dfile.all.1D'[2]' -stim_base 7
-stim_file 8  dfile.all.1D'[3]' -stim_base 8
-stim_file 9  dfile.all.1D'[4]' -stim_base 9
-stim_file 10 dfile.all.1D'[5]' -stim_base 10

-iresp 1 TMirf -iresp 2 HMirf
-iresp 3 TPirf -iresp 4 HPirf
-full_first -fout -tout -nobout -xjpeg Xmat
-bucket ${subj}_func

```

Movement  
regressors-of-  
no-interest:  
output from  
3dvolreg

- Output HRF (**-iresp**) 3D+time dataset for each stimulus class
  - Each of these 4 datasets will have TR=1.0 s and have 15 time points ( $\beta$  weights for lags 0..14)
  - Can plot these HRF datasets atop each other using **Dataset#N** plugin
  - Useful for visual inspection of regions that GLTs tell you have different responses for different classes of stimuli
- **-nobout** = don't output statistics of baseline parameters
- **-bucket** = save statistics into dataset with this prefix
- **-xjpeg** = save an image of the **R** matrix into file **Xmat.jpg**

## Script for Deconvolution - 9

```

-gltsym ../misc_files/contrast1.1D -glt_label 1 FullF
-gltsym ../misc_files/contrast2.1D -glt_label 2 HvsT
-gltsym ../misc_files/contrast3.1D -glt_label 3 MvsP
-gltsym ../misc_files/contrast4.1D -glt_label 4 HMvsHP
-gltsym ../misc_files/contrast5.1D -glt_label 5 TMvsTP
-gltsym ../misc_files/contrast6.1D -glt_label 6 HPvsTP
-gltsym ../misc_files/contrast7.1D -glt_label 7 HMvsTM

```

- Run many GLTs to contrast various pairs and quads of cases
- New feature: **-gltsym** = specify  $\beta$  weights to contrast using **-stim\_label** names given earlier on the command "line"
  - Simpler than counting 0s and  $\pm 1$ s to fill out GLT matrix numerically
- Example: file **contrast2.1D** is the single line below:
 

```
-ToolMovie +HumanMovie -ToolPoint +HumanPoint
```

 which means to put "-1" in the matrix for all 15 lags for stimuli #1 and #3 and "+1" in the matrix for all 15 lags for stimuli #2 and #4
  - This is the "Human vs Tools" contrast (labeled **HvsT** via **-glt\_label**)
  - Sum of the 30 "Tool"  $\beta$  weights subtracted from Sum of the 30 "Human"  $\beta$  weights
  - Testing: % signal change for Human stimuli different than Tool stimuli?

## Script for Deconvolution - 10

```

3dbucket -prefix ${subj}_func_slim -fbuc \
    ${subj}_func+orig'[0,125..151]'

foreach cond (TM HM TP HP)

    3dTstat -prefix ${subj}_${cond}_irf_mean \
        ${cond}irf+orig

    adwarp -apar ${subj}spgr+tlrc -dxyz 3 \
        -dpar ${subj}_${cond}_irf_mean+orig
end

cd ..
end

#=====
# End of script!
# Take ${subj}_${cond}_irf_mean+tlrc datasets into 3dANOVA3
#=====

```


Extract a subset of interesting statistics sub-bricks into a “slimmed-down” functional dataset

Compute HRF means across all lags 0..14 for each of the 4 stimuli types

Transform this individual's mean % signal results into Talairach coordinates for group analyses

End of loop over subjects; go back to upper directory whence we started

## Results: Humans vs. Tools



• Color overlay is **HvsT** contrast

• **Blue** (upper) curves: Human HRFs

• **Red** (lower) curves: Tool HRFs


## Yet More Fun 3dDeconvolve Options

- **-mask** = used to turn off processing for some voxels
  - ★ speed up the program by not processing non-brain voxels
- **-input1D** = used to process a single time series, rather than a dataset full of time series
  - ★ test out a stimulus timing sequence
  - ★ **-nodata** option can be used to check for collinearity
- **-censor** = used to turn off processing for some time points
  - ★ for time points that are “bad” (e.g., too much movement)
- **-sresp** = output standard deviation of HRF estimates
  - ★ can plot error bands around HRF in AFNI graph viewer
- **-errts** = output residuals (i.e., difference between fitted model and data)
  - ★ for statistical analysis of time series noise
- **-jobs *N*** = run with multiple CPUS — *N* of them
  - ★ extra speed, if you have a dual-CPU system (or more)!

## 3dDeconvolve with Free Timing

- The fixed-TR stick function approach doesn’t work well with arbitrary timing of stimuli
  - ★ When subject actions/reactions are self-initiated, timing of activations cannot be controlled
- If you want to do deconvolution (vs. fixed-shape analysis), then must adopt a different basis function expansion approach
  - ★ One that has a finite number of parameters but also allows for calculation of ***h(t)*** at any arbitrary point in time
- Simplest set of such functions are closely related to stick functions: **tent functions**


$$T(x) = \begin{cases} 1 - |x| & \text{for } -1 < x < 1 \\ 0 & \text{for } |x| > 1 \end{cases}$$



## Tent Functions = Linear Interpolation

- Expansion in a set of spaced-apart tent functions is the same as linear interpolation

$$\beta_0 \cdot T\left(\frac{t}{L}\right) + \beta_1 \cdot T\left(\frac{t-L}{L}\right) + \beta_2 \cdot T\left(\frac{t-2 \cdot L}{L}\right) + \beta_3 \cdot T\left(\frac{t-3 \cdot L}{L}\right) + \dots$$



- Tent function parameters are also easily interpreted as function values (e.g.,  $\beta_2$  = response at time  $t = 2 \cdot L$  after stim)
- User must decide on relationship of tent function grid spacing  $L$  and time grid spacing TR (usually would choose  $L \geq \text{TR}$ )
- Fancy name for tent functions: piecewise linear B-splines

## Tent Functions: Average Signal Change

- For input to group analysis, usually want to compute average signal change
  - ★ Over entire duration of HRF (usual)
  - ★ Over a sub-interval of the HRF duration (sometimes)
- In previous slide, with 6  $\beta$  weights, average signal change is

$$\frac{1}{2} \beta_0 + \beta_1 + \beta_2 + \beta_3 + \beta_4 + \frac{1}{2} \beta_5$$


- First and last  $\beta$  weights are scaled by half since they only affect half as much of the duration
- In practice, may want to use  $0 \cdot \beta_0$  since immediate post-stimulus response is not hemodynamically correct
- $\beta$  weights are output into the “bucket” dataset produced by **3dDeconvolve**
- Can then be combined into a single number using **3dcalc**

## 3dDeconvolve -stim\_times


- Direct input of stimulus timing, plus a response model
- Specifies stimuli, instead of using `-stim_file`
- `-stim_times k tname rtype`
  - ★ `k` = stimulus index (from 1 to `-num_stimts` value)
- `tname` = name of `.1D` file containing stimulus times (seconds)
  - ★ **N.B.:** TR stored in dataset header must be correct!
- `rtype` = name of response model to use for each stimulus time read from `tname` file
  - ★ **GAM** = gamma variate function from **waver** (fixed-shaped analysis)
  - ★ **TENT(b, c, n)** = tent function deconvolution, ranging from time `s+b` to `s+c` after each stimulus time `s`, with `n` basis functions (divided evenly over `c-b` seconds, into `n-1` intervals)
  - ★ several other `rtype` options available (experimental)
- Can mix `-stim_file` and `-stim_times` as needed
  - ★ e.g., movement parameter regressors at each TR

## Two Possible Formats of Timing File

- A single column of numbers
  - ★ One stimulus time per row
  - ★ Times are relative to first image in dataset being at `t=0`
  - ★ May not be simplest to use if multiple runs are catenated
- One row for each run within a catenated dataset
  - ★ Each time in  $j^{\text{th}}$  row is relative to start of run  $\#j$  being `t=0`
  - ★ If some run has NO stimuli in the given class, just put a single “\*” in that row as a filler
    - Different numbers of stim per run are OK
    - At least one row must have more than 1 time  
(so that this type of timing file can be told from the other)
- Two methods are available because of users’ diverse needs
  - ★ **N.B.:** if you chop first few images off the start of each run, the inputs to `-stim_times` must be adjusted accordingly



```
4.7
9.6
11.8
19.4
```



```
4.7 9.6 11.8 19.4
*
8.3 10.6
```




## Other Recent-ish Upgrades

- See <http://afni.nimh.nih.gov/doc/misc/3dDeconvolveSummer2004/>
- Equation solver: Gaussian elimination to compute **R** matrix pseudo-inverse was replaced by SVD (like principal components)
  - ★ Advantage: smaller sensitivity to computational errors
  - ★ “Condition number” and “inverse error” values are printed at program startup, as measures of accuracy of pseudo-inverse
  - ★ Condition number < 1000 is good
  - ★ Inverse error < 1.0e-10 is good
- **3dDeconvolve\_f** program can be used to compute in single precision (7 decimal places) rather than double precision (16)
  - ★ For better speed, but with lower numerical accuracy
  - ★ Best to do at least one run **both** ways to check if results differ significantly (SVD solver *should* be safe)

## Recent Upgrades - 2

- New **-xjpeg xxx.jpg** option will save a JPEG image file of the columns of the **R** matrix into file **xxx.jpg** (and an image of the pseudo-inverse of **R** into file **xxx\_psinv.jpg**)

Constant and linear baselines for each run (**-polort 1**)



Simple regression functions created by **waver** and input by **-stim\_file**

Why 'x' instead of 'R'? Because SPM calls this the 'X' matrix, not the 'R' matrix.

## Recent Upgrades - 3

- Matrix inputs for **-glt** option can now indicate lots of zero entries using a notation like **30@0 1 -1 0 0** to indicate that 30 zeros precede the rest of the input line
  - ★ Example: 10 imaging runs and **-polort 2** for baseline
  - ★ Can put comments into matrix and .1D files, using lines that start with **#** or **//**
  - ★ Can use **\** at end of line to specify continuation
- Matrix input for GLTs can also be expressed symbolically, using the names given with the **-stim\_label** options:
 

```
-stim_label 1 Ear -stim_maxlag 1 4
-stim_label 2 Wax -stim_maxlag 2 4
```

  - ★ Old style GLT might be
 

```
{zeros for baseline} 0 0 1 1 1 0 0 -1 -1 -1
```

Sum of Ear – Sum of Wax (lags 2..4)
  - ★ New style (via **-gltsym** option) is
 

```
Ear[2..4] -Wax[2..4]
```

## Recent Upgrades - 4

- New **-xsave** option saves the **R** matrix (and other info) into a file that can be used later with the **-xrestore** option to calculate some extra GLTs, without re-doing the entire analysis (goal: save some time by not recomputing)
- **-input** option now allows multiple 3D+time datasets to be specified to automatically concatenate individual runs into one file 'on the fly'
  - ★ Avoids having to use program **3dTcat**
  - ★ User must still supply full-length **.1D** files for the various input time series (e.g., **-stim\_file**, **-stim\_times**)
  - ★ **-concat** option will be ignored if this option is used
    - Break points between runs will be taken as the break points between the various **-input** datasets
- **-polort** option now uses Legendre polynomials instead of simple  $1, t, t^2, t^3, \dots$  basis functions (more numerical accuracy)

## Recent Upgrades - 5

- **3dDeconvolve** now checks for duplicate **-stim\_file** names and for duplicate matrix columns, and prints warnings
  - ★ These are not fatal errors
    - If the same regressor is given twice, each copy will only get half the amplitude (the “beta weight”) in the solution
- All-zero regressors are now allowed
  - ★ Will get zero weight in the solution
    - A warning message will be printed to the terminal
  - ★ Example: task where subject makes a choice for each stimulus (e.g., male or female face?)
    - You want to analyze correct and incorrect trials a separate cases
    - What if a subject makes no mistakes? Hmmm...

## Recent Upgrades - 6

- Recall: **-iresp** option outputs the HRF model for one stimulus
  - ★ When used with **-stim\_times**, values are usually output using the dataset TR time spacing
  - ★ Can changes to a different grid via new **-TR\_times dt** option, which sets the output grid spacing for **-iresp** to **dt** for HRF models computed via **-stim\_times**
    - Is useful for producing nice smooth pictures of HRF
    - Also works with **-sresp** option (= std.dev. of HRF)
- **Difficulty**: using GLTs with results from **-stim\_times**
  - ★ GLTs operate on regression coefficients
  - ★ For advanced (experimental) **rtype** models, regression coefficients don't correspond directly to HRF amplitudes
    - Exceptions: **GAM**, **TENT**, **BLOCK**


## Upgrades – *Planned* or *Dreamed of*

- Automatic baseline normalization of input time series
- Automatic mask generation (à la **3dAutomask** program)
- Spatial blur (à la **3dmerge -lblur**)
- Time shift input before analysis (à la **3dTshift** program)
- Negative lags for **-stim\_file** method of deconvolution
  - ★ for pre-stimulus cognition/anticipation
  - ★ **-stim\_times** already allows pre-stimulus response
- ‘Area under curve’ addition to **-gltsym** to allow testing of pieces of HRF models from **-stim\_times**
- Slice- and/or voxel-dependent regressors
  - ★ For physiological noise cancellation, etc.
- Floating point output format
  - ★ Currently is shorts + scale factor

## Advanced Topics in Regression

- Can have activations with multiple phases that are not always in the same time relationship to each other; e.g.:
    - a) subject gets cue #1
    - b) variable waiting time (“hold”)
    - c) subject gets cue #2, emits response
      - which depends on both cue #1 and #2
- } timing of events is known
- ★ Cannot treat this as one event with one HRF, since the different waiting times will result in different overlaps in separate responses from cue #1 and cue #2
  - ★ Solution is multiple HRFs: separate HRF (fixed shape or deconvolution) for cue #1 times and for cue #2 times
    - o Must have significant variability in inter-cue waiting times, or will get a nearly-collinear model
      - impossible to tell tail end of HRF #1 from the start of HRF #2, if always locked together in same temporal relationship
    - o How much variability is “significant”? Good question.

## Even More Complicated Case

- Solving a visually presented puzzle:
  - a) subject sees puzzle
  - b) subject cogitates a while
  - c) subject responds with solution

timing of events is measured
- The problem is that we expect some voxels to be significant in phase (b) as well as phases (a) and/or (c)
- Variable length of phase (b) means that shape for its response varies between trials
  - ★ Which is contrary to the whole idea of averaging trials together to get decent statistics (which is basically what linear regression amounts to, in a fancy sort of way)
- Could assume response **amplitude** in phase (b) is constant across trials, and response **duration** varies directly with time between phases (a) and (c)
  - ★ Need three HRFs; phase (b)'s is a little tricky to generate using waver, but it could be done

## Noise Issues

- “Noise” in FMRI is caused by several factors, not completely characterized
  - ★ MR thermal noise (well understood, unremovable)
  - ★ Cardiac and respiratory cycles (partly understood)
    - In principle, could measure these sources of noise separately and then try to regress them out
      - ➡ RETROICOR program underway (R Birn & M Smith of FIM/NIMH)
  - ★ Scanner fluctuations (e.g., thermal drift of hardware)
  - ★ Small subject head movements (10-100 mm)
  - ★ Very low frequency fluctuations (periods longer than 100 s)
- Data analysis should try to remove what can be removed and allow for the statistical effects of what can't be removed
  - ★ “Serial correlation” in the noise time series affects the  $t$ - and  $F$ -statistics calculated by **3dDeconvolve**
  - ★ At present, nothing is done to correct for this effect (by us)

## Nonlinear Regression

- Linear models aren't everything
  - ★ e.g., could try to fit HRF of the form  $h(t) = a \cdot t^b \cdot e^{-t/c}$
  - ★ Unknowns  $b$  and  $c$  appear nonlinearly in this formula
- Program **3dNLFim** can do nonlinear regression (including nonlinear deconvolution)
  - ★ User must provide a C function that computes the model time series, given a set of parameters (e.g.,  $a$ ,  $b$ ,  $c$ )
  - ★ Program then drives this C function repeatedly, searching for the set of parameters that best fit each voxel
  - ★ Has been used to fit pharmacological wash-in/wash-out models (difference of two exponentials) to fMRI data acquired during pharmacological challenges
    - e.g., injection of nicotine, cocaine, ethanol, etc.
    - these are tricky experiments, at best