

## AFNI Jazzercise

Please read the following questions and use your AFNI know-how to answer them. Hints to answering these questions are available in the “Hints” handout. The answers to these questions can be found in the “Answers” handout.

1. The dataset **AFNI\_data3/afni/func\_slim+orig** contains 7 sub-bricks of statistical data. Use **3dbucket** to create a smaller version of this dataset that contains only the sub-bricks: #0, 3-6. Name this new dataset **some\_stats**.  
*Why: To understand the layout of AFNI datasets and sub-bricks. You'll probably be interested only in specific sub-bricks and not all the sub-bricks that include baseline fit statistics.*
2. In directory **AFNI\_data3/afni** you will find three anatomical datasets: **anat1+orig**, **anat2+orig**, **anat3+orig**. These datasets are 3 separate anatomical scans of a single subject. They have already been aligned. Average them together into a single dataset called **anat\_mean+orig**. Notice that the result looks ‘cleaner’, since the noise has been reduced.  
*Why: averaging reduces the noise.*
3. Use AFNI’s two skull-stripping programs, **3dIntracranial** and **3dSkullStrip**, to remove the skull from dataset **AFNI\_data3/afni/anat+orig**. Name the output file from 3dIntracranial **anat\_3dIntra+orig** and the output file from 3dSkullStrip **anat\_3dSkull+orig**. Compare the two output datasets. Did one program do a better job at skull stripping or are the results similar? (Note: 3dSkullStrip may take a few minutes to run so be patient).  
*Why: Removing the skull is useful for image registration and creating a brain-specific mask.*
4. Creating and Playing with ROI Masks:
  - a. The dataset **AFNI\_data3/afni/func\_slim+orig** has beta values and t-stats for 2 stimulus classes, **fpos** and **fneg**. Use **3dcalc** to create a mask called **PN\_mask** that is **1** everywhere that both the fpos t-stat and the fneg t-stat values are greater than 4.2, and **0** everywhere else.  
*Why: Combining the results in a single mask is useful for a simple conjunction analysis.*
  - b. Similar to part a, create a conjunction mask that is **1** wherever  $a > 4.2$  (from fpos t-stat sub-brick), **2** wherever  $b > 4.2$  (from fneg t-stat sub-brick), **3** wherever both are true, and **0** otherwise. Name this dataset **PN\_mask\_4+orig** (since it contains 4 values).  
*Why: Using masks with values of powers of two are easier to look at all possible combinations in conjunction analysis.*
  - c. Use the afni GUI to display this mask, **PN\_mask\_4+orig**, so that each mask value gets its own color. What does each color mean?  
*Why: The various combinations will each show a unique color that is easily visible and understood.*
  - d. Use **3dROIstats** to store the average time series from **epi\_r1+orig** into the text file **PN\_mean.1D**, where the mean is over the voxels in the mask (from part a), **PN\_mask+orig**.  
*Why: Mean time series curves for each ROI are useful for further analysis and display with other AFNI programs or with external applications like Excel and Matlab.*

5. Fun with 1D files:

- a. Create three 1-column files with the numbers 1-10 in one column of the first file, 11-20 in the second file, and 21-30 in the third file. (note: you might use 2 different AFNI programs to create each file)
- b. Catenate these 3 files into one 3-column file. Call this 1D file **3\_cols.1D**.
- c. Create a new file that contains columns 1, 2, 3, 3, 2,1, from part b (i.e., there will be a total of 6 columns in this new 1D file). Call this new 1D file **6\_cols.1D**.
- d. Now take the 6 columns from question 7b and average them together to create a new file with a single column. Call that new file **ex\_mean.1D**.

*Why: AFNI 1D programs can help to combine and analyze data from multiple voxels or masks.*

6. Fun with the AFNI GUI:

- a. Open **AFNI\_data3/afni/anat+orig** and in any one of the views (sagittal, axial, or coronal), change the gray-scale intensity range to be 500 minimum and 1500 maximum.
- b. Open **AFNI\_data3/afni/func\_slim+orig** and set the Full-F as the OLay and Threshold. Set the Threshold to F=8.0. Show only Positive values and set the color scale to show only 8 colors. Edit the color scale so that F-values between 14 and 28 are shown in lime green.
- c. View the above settings you created from part b in a sagittal slice. Make a jpeg file from sagittal slice #166 and name it **cool\_slide**.
- d. Switch to Talairach view and Talairach to the **right fusiform gyrus**.
- e. Change the display to show 6 sagittal slices all at once, in a 3x2 montage.
- f. Can you find the AFNI Mission statement hidden in the AFNI GUI?

*Why: The AFNI GUI has the flexibility to look at your data, and it's fun once you learn how.*

7. Doing Calculations in AFNI:

- a. Determine what type of data (short, float, etc) makes up dataset **AFNI\_data3/afni/func\_slim+orig**.

*Why: It's important to know how to find information about your dataset.*

- b. Calculate  $22.3 * 44.5$  using the simple calculating program in AFNI.

*Why: AFNI provides calculation tools including simple ones that are useful in scripts.*

8. Image Filtering:

- a. Smooth **AFNI\_data3/afni/epi\_r1+orig** with a 8mm FWHM filter. Name the output file **ex\_blur8**.
- b. Enhance **AFNI\_data3/afni/anat+orig** by emphasizing the minimum-valued voxels across +/- 3 voxels in the sagittal (z) direction. Name the output dataset **ex\_minz3**.

- c. Enhance dataset **ex\_minz3+orig** from part b by removing the noise with program **3danisosmooth**. Name the output dataset **ex\_aniso**. Use the **-viewer** option in this program to select the number of noise-removing iterations.

*Why: Image filtering provides some powerful ways to enhance the data. Data can be changed in radically different ways. Use this power wisely.*

9. Random Exercises with AFNI Datasets:

- a. Open dataset **AFNI\_data3/afni/anat+orig** dataset and find the spatial storage order (i.e., xyz-orientation). Re-orient it to LPI orientation and name the new output dataset **exLPI**.

*Why: Data may be required by specific programs to match other data or to match a specific orientation.*

- b. Open dataset **AFNI\_data3/afni/func\_slim+orig** and create 2 separate datasets: one with the 3<sup>rd</sup> sub-brick only and one with the 4<sup>th</sup> sub-brick only. Call the former dataset **ex\_fneg\_coef** and the latter **ex\_fneg\_tstat**.

*Why: Extracting specific data is useful for exporting to other software or for making the data fit in memory more easily if it's a large dataset.*

- c. Combine **ex\_fneg\_coef+orig** and **ex\_fneg\_tstat+orig** from part b into a single dataset called **ex\_fneg**, having the fneg Coef for sub-brick 0, and fneg t-stat for sub-brick 1.

*Why: Datasets can be manipulated to include only what you're interested in.*

- d. Convert dataset **AFNI\_data3/afni/func\_slim+orig** to Talairach coordinates with a 4mm<sup>3</sup> resolution. Use dataset **anat+tlrc** in the same directory as the data parent to perform the transformation on **func\_slim+orig**. Name the output file **func\_slim4mm**.

*Why: Talairach data will be 1mm<sup>3</sup> by default. This resolution is often not necessary because it doesn't reflect the resolution of the EPI data. It makes processing slower and taxes memory too.*

- e. Locate the maximum "Full-F" stat voxel value in dataset **func\_slim4mm+tlrc** and find the name of the Talairach atlas region that corresponds to that voxel's position.

*Why: Finding maximum activation can be scripted or searched interactively. Atlas regions should be used as a guide. The AFNI GUI includes various atlases that cite the regions associated with any specific voxel.*

- f. Dataset **AFNI\_data3/afni/anat+orig** was acquired sagittally and contains 124 slices. Create a new dataset that contains only slices 40-90 of anat+orig. Provide the new dataset with the prefix name **anat\_40\_90**.

*Why: If memory or processing speed is a constraint, you can work with only part of the data.*