
# REAL-TIME FMRI: setup, image monitoring, statistics, and feedback

# Ziad S Saad, PhD SSCC / NIMH & NINDS / NIH / DHHS / USA / EARTH













# Why bother?

- Image quality control
  - Spikes, distortion, ghosting, noise, ...
  - Amount of motion

Cox, RW et al. 95, Cohen, MS et al. 98, Frank, J. et al 99, Voyvodic, J. 99

Functional localization

Weiskopf, N. et al 04

- Localizer prior to main FMRI experiment for BCI or high-res imaging
- Pre operative scanning
- As Q/A in clinical settings or difficult / rare subject population
- 'scan to criteria'

Yang, S. et al 08

Weiskopf, N et al. 2007

- Teaching
- Feedback and Biofeedback

Reduce motion

Alter/interfere brain function

Yang, S. et al. 05

deCharms, RC, et al. 04 deCharms. RC. et al. 05

Control of task/ stimulus computer Posse S. et al. 03

LaConte SM. et al. 07 Yoo S. et al. 04

Classification/BCI

Signals in vegetative state


Owen AM et al 06

# **Outline**

- This talk will focus on AFNI's interface for realtime FMRI
  - A brief intro to the interactive interface
  - Demo I: simple image monitoring
  - Demo II: Demo I + GLM
  - Demo III: Feedback
  - Demo IV: Classification ?
  - AFNI & SUMA Automation


- Image quality control
  - Spikes, distortion, ghosting, noise, ...
  - Amount of motion


Cox, RW et al. 95, Cohen, MS et al. 98, Frank, J. et al 99, Voyvodic, J. 99 Weiskopf, N et al. 2007



- Image quality control
  - Spikes, distortion, ghosting, noise, ...
  - Amount of motion


Cox, RW et al. 95, Cohen, MS et al. 98, Frank, J. et al 99, Voyvodic, J. 99 Weiskopf, N et al. 2007





Real-time
Estimation
of
Functional
Activation

Real-time
Estimation
of
subject
movement







Fig. 1. From Weiskopf, N. et al. MRI 07


# Reduce Motion with Feedback

Feedback and Biofeedback

Yang, S. et al. 08

Reduce motion





# Activation in Vegetative State

Patient and control responses to audio instructions






Fig.1 from Owen AM et al, Science 06

# Classification

- Classification maps high dimensional pattern into a set of classes
  - This allows a complex brain activation pattern to be identified with a set of classes or brain states.
  - Useful in to providing intuitive feedback from activation of multiple areas
  - Useful for inferring brain state
    - Strother, Cox Savoy, Haxby, Beauchamp, Kendrick Kay's



# Brain Computer Interface



# Why bother?

### Reviews:

- Weiskopf N et al.: Real-time functional magnetic resonance imaging: methods and applications. *Magnetic Resonance Imaging 25 (2007)*
- Yang S et al.: Real-Time Functional Magnetic Resonance Imaging and its Applications. in *Brain Mapping Research* Developments, Bakker LN ed., Nova Publishing, New Jersey (2008)
- deCharms RC: Applications of real-time fMRI. Nature Reviews
   Neuroscience 9 (2008)
- deCharms RC: Reading and controlling human brain activation using real-time functional magnetic resonance imaging. *Trends* in Cognitive Sciences 11 (2007)

# The players

**Image Monitor** 


Scanner

AFNI

Real Time Setup

Stimulus Display

Real Time Receiver



### Real Time Setup

- A user-supplied set of commands that tell AFNI what to do with incoming data
- Can be done from shell commands or from within C code
- Communicates with AFNI through TCP/IP socket
- Sets up ROIs for AFNI\*. (need to show example of ROI selection.
   Either single subject-based or atlas based).

# Real time setup example

- A module from the demo
- @fast\_roi

```
@fast_roi -region CA_N27_ML::Hip \
    -region CA_N27_ML::Amygda \
    -base TT_N27_r2+tlrc. \
    -anat doe_SurfVol_Alnd_Exp+orig. \
    -roi_grid blur_vr_run1_motor_AFB003+orig. \
    -prefix hip_amy -time
```

Freesurfer based selection

# ROI selection options


- Standard atlases
  - TT\_Daemon :
    - Created by tracing Talairach and Tournoux brain illustrations.
    - Contibuted by Jack Lancaster and Peter Fox of RIC UTHSCSA
  - CA\_N27\_MPM, CA\_N27\_ML, CA\_N27\_PM :
    - Anatomy Toolbox's atlases, some created from cytoarchitectonic
    - studies of 10 human post-mortem brains
    - contributed by Simon Eickhoff, Katrin Amunts and Karl Zilles of IME,
       Julich,
- FreeSurfer, subject-based
- Functional localizer
- Etc.

## Standard-space atlas\_ROI selection

```
@fast_roi -region CA_N27_ML::Hip \
    -region CA_N27_ML::Amygda \
    -base TT_N27_r2+tlrc. \
    -anat doe_SurfVol_Alnd_Exp+orig. \
    -roi_grid blur_vr_run1_motor_AFB003+orig. \
    -prefix hip_amy -time
```

- less than 1min including skull stripping and transform to TLRC
- A couple of seconds for generating more ROIs

## Atlas-based ROIs




1- Strip skull2- Find xform to atlas space (about 40 secs, 2.5GhZ cpu)



3- Identify ROIs
4- Xform ROIs to native space
(about 2 seconds)



## Subject-based Anatomical ROIs



# The players Scanner Real Time Setup AFNI Real Time Receiver Stimulus Display

- Scanner
  - A user-supplied machine to acquire and reconstruct images in real time


# Scanner Real Time Setup Image Monitor AFNI Stimulus Display Real Time Receiver

- Image Monitor
  - An AFNI- or user- supplied program to wait for scanner reconstructed images
    - AFNI-supplied programs monitor files only:
      - Imon (Monitors GE's old dreaded I files)
      - Dimon (Monitors GE's DICOM images)
      - RTfeedme (Breaks up and AFNI readable dataset and sends it to AFNI)
    - User-supplied programs usually interface with scanner software directly
      - SIEMENS TRIO(?) via functors (S. LaConte BCM, E. Stein NIDA)
    - Usually only program that runs on scanner computer
      - But even that is not necessary
  - Image Monitor sends images or volumes to AFNI over TCP/IP socket as they become available

    Z.S.S 8-09

# Image Monitor example

- A module from the demo
- If there is time, add flowchart to point to small but pesky problems that they should watch for.



### AFNI

- AFNI-GUI application expects incoming images/volumes and processes them per the setup instructions
  - Assemble images/volumes into time series
  - Perform image registration
  - Perform (multi\*) linear regression
  - Send results to Real Time Receiver through TCP/IP socket
    - Raw, volume registered, or residual volume\*
  - Send raw or processed volumes to plugin that register to receive them
    - Much faster than TCP/IP (just a data pointer is passed)
    - Users typically write their own plugins (S. LaConte 3dSVM)
    - Plugins can also communicate with Real Time Receiver

Sand data from are defined POIs to Peal Time Possiver

# AFNI example

- Should probably show the kinds of data to AFNI sends
- Maybe a picture of a plugin
- Maybe some code to show how plugin requests data?
- Two tables of ROI-based values as they get sent

# Scanner Real Time Setup Image Monitor AFNI Stimulus Display Real Time Receiver

### Real Time Receiver

- AFNI- or User- supplied application that expects incoming data from AFNI and acts on it
  - Motion parameters
  - ROI-based data, all values or just average
  - Entire volumes of raw, or preprocessed data
  - Data from any RT plugin such as 3dsvm
- Process incoming data to your liking
- Optionally forward results to Stimulus Display either by serial connection, or TCP/IP\*

# Receiver example

Example from demo

- Speak of importance of scanner independent testing platforms
- Extend hand for collaboration

# "Help" sources

- Which programs are available for help
- Readme files

# Strategy for Manipulating Activation

Adapted from deCharms RC. TCS 07


- Providing strategy may be critical
  - Subjects overestimate ability to control activation
  - Start by providing strategy that activates ROIs regions providing feedback
- See literature on control of various areas
  - Somatomotor cortex

From LaConte S. – FMRI Advanced Issues ISMRM 09

- Posse 2001, Yoo 2002, deCharms 2004, Yoo 2004
- Parahippocampal place area
  - Weiskopf 2004
- Amygdala
  - Posse 2003
- Insular cortex
  - Caria 2007
- Anterior cingulate cortex
  - Weiskopf 2003, Yoo 2004, Birbaumer 2007, deCharms 2005

If incidental to task, minimize interference

Too much information!



• If incidental to task, minimize interference

Minimum Task Interference

**Enough information** 




Fig.7 from Yang, S. et al. Neuroimage 05

Fig.3 from Yang, S. et al. Neuroimage 05

- Make it appealing to subject
  - Turns out few get excited about graphs!
  - Fire on the beach = much more exciting



Figure 1d from deCharms RC. Nature 08

- OMG! Asteroids!
  - Keeps subject interested



History trace helps subject cope with FMRI response lag






# Interface Design







# Interface Design



From S. LaConte

ISMRM 09


# Interface Design



From S. LaConte

ISMRM 09

# Interface Design



From S. LaConte

ISMRM 09

# What to feedback?

- Which signal to use?
  - From original time series
  - From filtered\* time series
  - From regression (Beta/T/R) analysis
- Typically from one or more ROIs
  - Anatomical Atlas based
  - Single subject anatomy based
  - Group function based
  - Single subject localizer
- Combining information from multiple ROIs
  - Encode signals in VR scene

deCharms RC. 08

- Classifiers (ROI or whole brain), if models are known LaConte SM. 07
- What about noise confounds?
  - Control for respiration/cardiac with real-time RETROICOR\*
  - Include other physiological covariates in real-time\*
  - Include real-time baseline modeling

# <u>Acknowledgments</u>

Robert Cox Rick Reynolds

Stephen LaConte Tomas Ross

Julien Doyon