

Musculoskeletal MRI: Practical Protocols

Timothy J. Mosher, MD
Chief, Musculoskeletal Imaging and MRI
Penn State University College of Medicine
Hershey, PA



Learning Objectives

- Discuss general considerations in designing clinical MRI protocols
- Understand the role of contrast resolution in the context of imaging connective tissues
- Present clinical MRI protocols for MSK MRI at 3.0 T

What is a “Practical” MRI Protocol


- A practical protocol is not a perfect protocol
- The final product must satisfy many different customers with competing interests
- For every protocol there is a colleague with a better protocol


Organizational Guidelines

- Use detailed MRI requests and patient questionnaires to extract sufficient history
 - Mark site of tenderness/mass with fiducial marker
- Use targeted MRI protocols
 - Primary objective
 - Secondary objective
- Limit patient examination times to 45 minutes or less
- Always perform the most important sequence first
- Invest in education for technologists

General Approach for Designing Clinical MRI Protocols

1. Obtain sufficient signal to noise (SNR) to get the job done
2. Optimize contrast for the tissue that you are evaluating
3. Select image plane and resolution based on the anatomy that you are evaluating
4. Adjust acquisition parameters to minimize artifact






Principle 1: *Musculoskeletal MR imaging protocols are optimized for contrast resolution at the expense of imaging speed*

Optimizing Contrast Resolution




Pulse Sequence Selection



Commonly used sequences for MSK MR imaging


- Conventional spin echo
- Fast (Turbo) spin echo
- Gradient echo

Conventional Spin Echo




Conventional Spin Echo

- Advantages
 - Excellent contrast
 - Experience
 - Validation studies
- Disadvantages
 - Long imaging times




3.0 T Proton density weighted spin echo

Fast (Turbo) Spin Echo




Fast (Turbo) Spin Echo




Fast (Turbo) Spin Echo

- Advantages
 - Efficient signal acquisition
 - Fast imaging times
 - Higher spatial resolution
 - Excellent contrast
- Disadvantages
 - Image blurring
 - Magnetization transfer
 - High SAR




1.5 T Turbo Spin Echo PD-weighted

Gradient Echo



Gradient Echo


- Advantages
 - Fast imaging times
 - 3D acquisitions are feasible
- Disadvantages
 - Moderate image contrast
 - Prone to artifact
 - Metal artifact



1.5 T Water Excited T1-weighted GRE

Principle 2: *First optimize contrast based on tissue type, then adjust resolution based on anatomy*


MRI Contrast in MSK Imaging



Effect of collagen on tissue contrast

- Efficient spin-spin (T2) relaxation
 - Tissue T2 is inversely related to collagen concentration
 - Tissue T1 is less dependent on collagen concentration
 - Anisotropic arrangement of collagen fibrils produces an orientation dependence of T2 (Magic angle effect)
- Magnetization transfer
 - Collagen is the dominant macromolecular component for magnetization transfer

Field Dependence of Relaxation Times



Duewell SH. *et al. Radiology* 1995; 196:551-555
Fullerton GD *et al Radiology* 1985; 155:433-435 (tendon)

Sequence Selection for MSK MRI


- | | |
|--|--|
| <p>Soft tissues</p> <ul style="list-style-type: none"> – Muscle/Fat/Bone Marrow <ul style="list-style-type: none"> • T1 FSE • T2 FSE with fat suppression • STIR | <p>Connective Tissues</p> <ul style="list-style-type: none"> – Menisci/Ligaments/Tendons <ul style="list-style-type: none"> • T1 SE • T1 or PD FSE • PD FSE with fat suppression – Articular Cartilage <ul style="list-style-type: none"> • PD FSE with or without fat suppression • Fat suppressed T1 spoiled gradient echo |
|--|--|

Principle 3: *In evaluation of connective tissue pathology tissue contrast will primarily be determined by:*

Pulse Sequence


Echo Time

How do tissue properties influence MRI contrast?




Mow VC, Procter CS, Kelly MA. *Biomechanics of articular cartilage*. In: Basic biomechanics of the locomotor system, Nordin M, Frankel VH eds., pp 31-58 Philadelphia Lea and Febiger.

Effect of collagen on cartilage T2




Inverse Correlation of Cartilage T2 and Polarized Light Microscopy




Courtesy of MT Nieminen, Beth Israel Deaconess Medical Center, Boston, MA


The Magic Angle Effect in evaluation of connective tissue




Dependence of Cartilage T2 on Collagen Fibril Orientation



Magic Angle Effect: Tendons



Magic Angle Phenomenon




Tailoring the MRI protocol for evaluation of connective tissues


- Most clinical MSK MRI requests are for evaluation of connective tissue pathology
- Primary indications
 - Shoulder: rotator cuff tear
 - Knee: meniscal tear
 - Ankle: tendon or ligament tear

MSK MRI protocol must be designed to accurately characterize connective tissue pathology, with the critical factor being identification of a surgical lesion (i.e. tear)


Effect of Connective Tissue Pathology on T2




T2 changes with tendon pathology




Choice of TE in Evaluation of Tendons and Ligaments



Evaluation of the rotator cuff: Short TE Sequence



Evaluation of the rotator cuff: *Long TE Sequence*



Normal


Tendinosis

Full
Thickness
Tear


Principle 4: *In evaluation of connective tissues two TE values are often needed:*

Short TE: high sensitivity, low specificity



Long TE: low sensitivity, high specificity



Magnetization Transfer:
A Probe for Collagen Content





Physics of Magnetization Transfer



Magnetization transfer decreases signal intensity of bound water

Magnetization Transfer With FSE Pulse Sequences



Yao, L, Gentili A, Thomas A. *Incidental Magnetization Transfer Contrast in Fast Spin-Echo Imaging of Cartilage. JMIR*, 6(1),180-184(1996).

Clinical Example: *Effect of Magnetization Transfer on visualization of cartilage*




T1-weighted fat suppressed GRE




PD-weighted fat suppressed FSE

Clinical Application of 3D T1-GRE Cartilage Imaging




1.5 T 3D Fat Suppressed T1-weighted GRE



1.5 T Water excited 3D T1-weighted GRE

3.0 T Wrist Imaging




3D Water Excited GRE




Fat Suppressed Proton Density TSE

Contrast versus resolution in visualization of superficial fibrillation




Fat-suppressed 3D FLASH
1mm, 512² matrix
0.1 mm³ voxel size




3D DESS
2 mm, 256² matrix
1.0 mm³ voxel size

Improved visualization of superficial cartilage lesion with FSE




Fat sat T1-weighted GRE




Fat sat Proton Density FSE


Tissue Contrast at the Articular Surface




Fibrillated Cartilage



MT improves visualization of the articular surface of cartilage




Grade I: Blistering




MRI Findings: Focal elevation in Cartilage T2 with or without superficial fibrillation

Diagnostic impact of MT depends on tissue type and location




Gradient Echo versus FSE



T1-weighted GRE

T1-weighted FSE


Gradient Echo versus FSE



T1-weighted GRE

T1-weighted FSE

Conventional SE versus FSE in evaluation of meniscal tear




Conventional SE Sequence
TR: 2000 ms
TE 20 ms

TSE Sequence (ETL: 3)
TR: 2000 ms
TE 12 ms

Why are meniscal tears less conspicuous on FSE?

- Magnetization Transfer decreases signal intensity of fluid within tear
- Blurring due to T2 modulation of the point spread function

FSE Blurring




Tips to minimize FSE blurring

- Reduce the echo train length (< 6)
- Reduce the time interval between echoes (inter-echo spacing)
- Increase spatial resolution
- Less effect on FSE T2 weighted images than on short TE images

Dynamic Range




Use of fat suppression to increase dynamic range of tissue contrast



FSE Sagittal PD

Sagittal PD with fat suppression

28 Year old female with prior ACL repair and knee pain



TR/TE: 2000 ms/15 ms

TR/TE: 2000 ms/36 ms


3.0 T

Fat Suppression

Summary of tissue contrast considerations

- Signal intensity changes in connective tissues are dominated by T2 effects of collagen on water
- Magnetization transfer is a critical mechanism of contrast at tissue interfaces
- Short TE is needed to detect changes in the collagen matrix (tendon degeneration)
- Long TE is needed to characterize free fluid (*diagnose tendon tear*)
- Rapid T2 decay results in image blurring with fast (turbo) spin echo sequences
- Fat suppression is useful to increase dynamic range

Composing practical MRI MSK protocols



The Shoulder Protocol


- Primary indication: Rotator cuff evaluation
- Secondary indications
 - Instability
 - Labrum
 - Capsule
 - Cartilage
 - Muscles
 - Marrow
 - Periarticular soft tissues

Optimizing TE: *The Shoulder MRI Protocol*

| Sequence | TR (ms) | TE(eff) (ms) | ST (mm) | Matrix | FOV (cm) |
|--------------------------------|---------|--------------|---------|-----------|----------|
| Axial FSE PD with FS (ET 6) | 4000 | 30 | 3 | 512 x 512 | 18 |
| Coronal Oblique FSE T1 (ET 3) | 500 | 15 | 4 | 512 x 512 | 18 |
| Coronal Oblique FSE T2 (ET 12) | 4000 | 60 | 4 | 512 x 512 | 18 |
| Sagittal Oblique FSE PD | 4000 | 100 | 4 | 512 x 512 | 18 |

3.0 T Protocol

43 year old male with abduction weakness



The Knee Protocol


- Primary indication: detection and characterization of meniscal tear
- Secondary indications
 - Ligamentous injury
 - Osteochondral pathology
 - Soft tissue inflammation

3.0T Knee MRI Protocol

| Sequence | TR (ms) | TE(eff) (ms) | ST (mm) | Matrix | FOV (cm) |
|-------------------------------------|---------|--------------|---------|-----------|----------|
| Fat Sat Axial FSE PD with DE (ET 5) | 4200 | 30 | 3 | 512 x 512 | 16 |
| Sagittal FSE PD (ET 5) | 2500 | 15 | 4 | 512 x 512 | 16 |
| Fat Sat Sagittal FSE T2 (ET 6) | 5500 | 45 | 4 | 512 x 512 | 16 |
| Fat Sat Coronal FSE PD (ET 5) | 4200 | 30 | 4 | 512 x 512 | 18 |

3.0 T Protocol

23 year old with post-traumatic ACL insufficiency




The Hip Protocol

- Primary indication: pain unresponsive to conservative management
 - AVN
 - Labral tear
 - Acetabular femoral impingement
 - Greater trochanteric bursitis

Problem: Need for large region of coverage with high spatial resolution

29 year old professional hockey player with chronic groin pain




3.0T Hip MRI Protocol

| Sequence | TR (ms) | TE(eff) (ms) | ST (mm) | Matrix | FOV (cm) |
|--------------------------------------|---------|--------------|---------|----------------------------|----------|
| Coronal FSE T1 (ET 3) | 800 | 12 | 5 | 512 x 512 | 34 |
| Coronal FSE T2 with FS (ET 16) | 5500 | 80 | 5 | 1024 x 1024 SENSE factor 2 | 34 |
| Coronal axial and sagittal PD (ET 5) | 2000 | 30 | 3 | 512 x 512 | 16 |
| Coronal PD with fat sat (ETL: 9) | 2200 | 40 | 3 | 512 x 512 | 16 |

Phased array body coil both hips

Paired surface coil symptomatic hip


55 year-old female with left hip and buttock pain with exercise, suspect piriformis syndrome



1.5 T Coronal SSF T2 phased array body coil

3.0 T Coronal FSE T2

3.0 T Hip MRI



Coronal PD

Coronal PD with fat suppression

24 year old female soccer player with chronic hip pain



MR Arthrogram requested for suspected labral tear

Take Home Points

- Trade-off of contrast, resolution, and speed
 - MSK protocols optimized for contrast resolution
- Optimize contrast first then adjust anatomy
 - Contrast is tailored for the tissue type
 - Resolution is tailored for the anatomy
- Connective tissue contrast is strongly influenced by collagen
 - Short T2
 - Magnetization transfer
 - Orientation dependence of signal intensity
- Dynamic range is set by fat, use of fat suppression allows shorter TE's to be used to obtain fluid sensitive sequences

