

MRA at 3.0T

Paul Finn, MD

David Geffen School of Medicine at UCLA



Acknowledgements

- Stefan Ruehm, MD
- Kambiz Nael, MD
- Roya Saleh, MD
- Anthony Ton, MD
- Mayil Krishnam, MD
- Sergio Godines, RT
- Glen Nyborg, RT
- Howard Dinh, MD
- Carissa Fonseca, Ph
- Gerhard Laub, PhD
- Vibhas Deshpande, PhD
- Hanns Weinmann, PhD

General Trends at 3.0T vs 1.5T

Increased SAR

- Limits flip angles and minimum TR for high-performance sequences
 - SSFP cine
 - Spin echo train imaging
 - CEMRA

B1 Inhomogeneity

- Shading in some regions
- Inhomogeneous contrast in some regions
- May make calibration of RF transmitter voltages difficult – varying flip angles within body regions

Carotid and Thoracic MRA at 3.0T


Carotid Disease

Bifurcations

Origins

Intracranial

Vertebro-basilar



CEMRA Procedure: typical

Center on head and do pre-Gd brain sequences

Center on neck via remote table movement (typical offset 140 mm)

Multiplanar scout images

timing-localizer (< 60 secs)

Run pre- and post-contrast, breath-held 3D MRA

Move table and do post-Gd brain sequences

Contrast Injection Scheme

Electronic infusion pump

Timing run with 2 ml Gd

MRA: 20 ml Gd at ~ 1.2 ml /sec flow rate


30 ml saline flush @ 1.2ml/sec

Also works at 1.5T!


Carr JC, Ma J, Deshpande V, Pereles FS, Laub GL, Finn JP. High-Resolution Breath-Hold Contrast-Enhanced MR Angiography of the Entire Carotid Circulation. *AJR* 2002; 178:543-549


Yang CW, MD; Carr JC, Futterer SF, Morasch MD, Benson P, Yang BP, Shors SM, Finn JP. Contrast-Enhanced Magnetic Resonance Angiography of the Vertebrobasilar Circulation. *AJNR* 2005 Sep;26(8):2095-101

TEST BOLUS



Timing-run /Localizer





Time-resolved 3D

- Nael K, Michael H, Villablanca VP, et al. Time-resolved 3D MRA of the Head and Neck at 3.0T: Initial Results. *Investigative Radiology* 2006;41:116-124.

UCLA

CEMRA Protocol: Time-resolved 3D

Fast 3D GRE sequence:

TR/TE 2.1/0.9 ms, FA 15°, BW 1200 Hz/pixel

Matrix 384 x 256

Partition thickness 4 mm


GRAPPA x 3

5 ml Gd at 3 ml/s followed by 30 ml of saline flush


Temporal resolution ~ 1.5 secs

UCLA

Time-Resolved 3D MRA



5ml Gd, 1.5 secs per 3D phase



UCLA

CEMRA Protocol: High Resolution 3D

Nael K, Villablanca P, Pope W, Laub G, Finn JP.
3D Contrast Enhanced MRA of the Head & Neck at 3.0T: Using Highly Acceleration Parallel Acquisition. *Radiology*, in press.

UCLA

CEMRA Protocol: High Resolution 3D

Fast 3D GRE sequence:

TR/TE 3/1.2 ms, FA 20°, BW 790 Hz/pixel

Matrix 460 x 576 (0.84 x 0.67 mm²)







Partition thickness 0.75 mm







GRAPPA x4

22 ml Gd at 1.2 ml/s followed by 30 ml of saline flush

A coronal 3D slab included the aortic arch, carotids and vertebro-basilar circulation in a 20 second breath-hold


UCLA





Perfusion Semi-quantitative Analysis

Gamma variate fit on a pixel by pixel basis calculates: MUS, MTT, TTP, MSI



UCLA

Conclusion


Anatomical vascular structures were seen with higher SNR at 3.0 T (consistent with theoretical predictions)

Parenchymal enhancement (often related to microvascular perfusion) is better seen at 1.5 T than 3.0 T


No significant difference was found between F.Q. or quantitative perfusion analysis at 1.5 & 3.0 T

UCLA


3.0T Dissection



3.0T Dissection: second phase




TIM Trio: - dissection




6ml Gd, 12 measurements each 1.7 s apart
21 s breath hold: iPAT x 3; TREAT


High Resolution MRA at 3.0 Tesla




High Resolution MRA at 3.0 Tesla




Renal Transplant: 3.0T Tim Trio



GI symptoms: M 55yrs



Bilateral Renal Aneurysms




Breath-hold, ~20 sec
Parallel imaging iPAT x 3
injection of 25 ml Gd-DTPA


TIM Trio: Branch Renal Artery Stenosis




TIM Trio: Renal Arteries




TIM Trio:— whole body



TIM Trio:— whole body



TIM Trio:— whole body



Conclusion

MRA at 3.0T is versatile and mature

Only modest doses of contrast are required

Attention to detail is important

Post-processing can greatly enhance visualization

Thank you

