

# Clinical MRI: From Physical Principles to Practical Protocols

## Cardiovascular Methods

Christine H. Lorenz, PhD  
Siemens Corporate Research  
Baltimore, MD

SIEMENS

2

## Disclaimers/Disclosures

1. The use of contrast media has not been approved by the FDA for use in Cardiovascular MR. Off label use of contrast media will be discussed in this presentation.
2. The information in this presentation is strictly educational and is not intended to be used for instruction as to the practice of medicine.
3. The speaker is an employee of Siemens Corporate Research

SIEMENS

3

## Outline

- Physiology and cardiac MR acquisitions
- Morphology
- Function
- Flow
- Perfusion
- Late enhancement

SIEMENS

4

## Outline

- Physiology and cardiac MR acquisitions
- Morphology
- Function
- Flow
- Perfusion
- Late enhancement

SIEMENS


5

## Motion in Cardiovascular MR

### Cardiac



### Respiratory





SIEMENS


6

## Motion Compensation Methods

- ECG Triggering
  - synchronize data acquisition over multiple heartbeats
- Respiratory gating
  - restrict acquisition to quiet periods of cycle
  - breath-holding
- Real time imaging
  - freeze motion



- SIEMENS 9
- ## MRI-Related ECG Noise Sources
- $B_0$  related (main static field)
    - magnetohydrodynamic effect
    - flux change due to motion (respiratory and cable)
  - $B_1$  related (radio-frequency field)
    - RF pulses
    - electrode and ECG amplifier effects
  - G related (gradient switching)



- SIEMENS 11
- ## Diagnostic Value of ECG in MR Environment
- Not Diagnostic
    - P-wave, ST-segment, T-wave, T-axis, J-point
    - PQ interval, QT-interval, P duration
  - Diagnostic
    - heart rate (accurate QRS detection required, or manual evaluation)
    - electrical axis of the heart
  - Heart rate changes due to  $B_0$
- Fischer SE et al ISMRM 2000

- SIEMENS 12
- ## Skin Preparation
- Remove chest hair if needed
  - Clean and roughen skin surface
    - Remove oil and perspiration to allow better contact of electrodes
    - Clean with abrasive prep pad or gel (NuPrep)
    - Dry area completely!!!

SIEMENS
13

## Use MR Compatible ECG Equipment

- MR compatible ECG system
  - Bruker
  - In Vivo
  - Medrad
- MR compatible ECG pads
  - check with MR vendor for recommendation
  - check expiration date
  - do not leave pads exposed to air-will dry out (baseline noise)
- Non MR pads
  - worst case ~ burns
  - RF interference ~ unable to scan!

SIEMENS
14

## Lead Placement Tips: Standard ECG

- The R wave with the highest amplitude will be attained by placing the positive and negative electrodes, **parallel to the electrical axis** of the heart
- The reference electrode can be placed in any position

SIEMENS
15

## Principle of the VCG (Vectorcardiogram)

- Two-channel ECG system
- Signals from two pairs of electrodes plotted against each other to remove time
- Resulting graph is spatial representation of electrical activity in the heart.

SIEMENS
16

## VCG Triggering Principle

Gradient noise, flow artifacts

Magnitude Signal

Acceptance Range

- R-wave is spatially separate from artifacts.
- 2D direction of R-waves is determined during learning-phase, while patient is outside magnet bore.
- Direction information is used to differentiate R-waves from artifacts.

US Patent 5,987,348 SE Fischer, SA Wickline and CH Lorenz

SIEMENS
17

## 2D VCG Lead Setup Philips Medical Systems

SIEMENS
18

## Wireless VCG System Siemens Medical Solutions

Fixed lead placement  
Single reference electrode  
One vertical lead (AVF), and one horizontal lead (Lead I)

SIEMENS19

## Summary: ECG Triggering

- Don't start scanning unless ECG is good
- Skin preparation is important for signal quality
- For non-VCG systems, placement of electrodes based on patient characteristics can greatly influence quality of ECG signal

SIEMENS20

## Electromechanical Activity in the Heart

ECG

Pressure (Aortic, LV, LA)

LV Volume

SIEMENS21

## Filling k-space Using ECG Triggering

ECG

R Wave Detection

Trigger Delay (TD)

k-Space

SIEMENS22

## ECG Triggered Acquisition with Respiratory Gating

k-space filled when respiratory motion acceptable and at correct point in cardiac cycle

SIEMENS23

## Navigator Gated Imaging Sequences

Short acquisition window in diastasis (<100 ms)

Navigator efficiency 30-50 % heartbeats accepted

Navigator pulse

Fat sat pulse or other contrast enhancing pulse

k-space

SIEMENS24

## Navigator Technique: Prospective Mode







- Detect respiratory motion
- Accept data only in expiration
- Compensate by adjusting slice position (slice following)

Philips Medical Systems

Siemens Medical Solutions

SIEMENS25

Outline

 Physiology and cardiac MR acquisitions
  Morphology
  Function
  Flow
  Perfusion
  Late enhancement

SIEMENS26

Imaging Heart Anatomy

- Clinical goals**
  - Determine relationship between various structures in the heart and in the thorax
  - Identify and characterize tumors, thrombus, fatty infiltration of myocardium, etc
  - Identify and characterize vessel wall abnormalities
- Sequence requirements**
  - Good contrast between blood and myocardium
  - Ability to have T1 or T2 weighting to characterize fat, tumors etc.
  - Moderate spatial resolution 1 x 1 mm in plane, 5-7 mm through plane
  - Good suppression of artifacts

SIEMENS27

Examples of Anatomic Imaging







SIEMENS28

Spin Echo Based Black-blood Techniques

	Spin echo	Turbo spin echo	Haste
Lines acquired per beat	1	8-32	64-128
Temporal resolution	10 – 40 ms	30 –300 ms	200 – 400 ms
Scan time	128 – 384 heartbeats	4 – 32 heartbeats	1 heartbeat

SIEMENS29

Breath-hold Black-blood Turbo Spin Echo




- Breath-hold to reduce respiratory artifacts.
- Diastolic gating to reduce cardiac motion artifacts.
- Double-inversion blood effectively eliminates signal from flowing blood.

Matrix 120x256, FOV 260 x 350mm  
6mm slice, 16 heart beats







SIEMENS30







Black-blood Turbo Spin Echo



with
 without
 with
 without

Blood signal artifacts in turbo spin echo can be severe without blood nulling







SIEMENS43

## Functional Imaging In The Heart

- Clinical goal
  - View and quantitate left and right ventricular function, wall motion, and wall thickening
- Sequence Requirements
  - Moderate spatial resolution (1-2 mm in-plane, 6-8 mm through plane)
  - High temporal resolution (< 50 ms/frame)
  - Suppression of respiratory motion (breathhold, real-time, respiratory gating)
  - High contrast between blood and myocardium

SIEMENS44

## Examples: Functional imaging in the Heart





Standard BH acquisition

Real-time acquisition

SIEMENS45

## Multiple Phase ECG Triggered Imaging (Cine Imaging)



Example: 3 phases


SIEMENS46

## Strategies for Shortening Scan Time

- Collect more than one k-space line per phase per beat (segmentation)
- Share collected k-space lines between phases (segmentation with echo sharing)
- Collect fewer lines than in full matrix and use properties of k-space or sensitivity maps to 'fill in' the rest
  - partial Fourier
  - parallel imaging (SENSE, SMASH, GRAPPA)
  - undersampled radial imaging
- Go really fast...get all lines at once (real time imaging)

SIEMENS47

## Improvement in Scan Time Cardiovascular MR









functional imaging: time for one cine loop

SIEMENS48

## Penalties for Shorter Scan Times

- Less spatial resolution
- Less temporal resolution
- Artifacts





SIEMENS 55

## Phase Contrast Cine Images

rephased

magnitude

phase

magnitude of flow compensated signal S1

magnitude of difference signal S2 – S1

phase of difference signal  $\gamma$


flow bright background visible

flow bright background suppressed


forward flow bright reverse flow black background midgray

SIEMENS 56

## Through-Plane vs In-Plane Encoding of Velocity



In-Plane Sagittal Aorta



Thru-Plane Axial Aorta

SIEMENS 57


## Typical Peak Velocities (cm/sec)

Pulmonary Artery	70-130
Aorta	100 – 175
Carotid Artery	80 – 120
External Iliac Artery	81 – 120
Carotid Syphon	55
Common Femoral Artery	115
Basilar Artery	40
Superficial Femoral Artery	90
Vertebral Artery	40
Popliteal Artery	70
Sagittal Sinus Vein	10
Peripheral Veins	5 – 10

SIEMENS 58


## What Happens If The VENC Is Wrong?

- VENC too low
  - aliasing
  - can correct to some extent, recommended to re-acquire with appropriate higher venc
- VENC too high
  - loss of 'velocity' resolution (poor dynamic range)




SIEMENS 59


## Effect of VENC Selection




VENC optimal  
(max black/white & no aliasing)



VENC only slightly too low  
(very little aliasing)




VENC much too high  
(washed-out black/white)




VENC much too low  
(very much aliasing)

SIEMENS 60


## Outline




Physiology and cardiac MR acquisitions




Morphology




Function



Flow



Perfusion



Late enhancement

SIEMENS61

# Myocardial Perfusion Imaging

- Clinical goals
  - Detect deficits in blood flow to the heart muscle at rest or under stress
- Sequence requirements
  - At least 1 image per heartbeat per slice
  - Slices to cover relevant segments of the heart (3-5)
  - High T1 sensitivity
  - Good SNR
  - Good spatial resolution (at least 2 x 2 in plane)
  - Resistance to motion artifacts and variations in heart cycle length

SIEMENS62

# Contrast Enhanced MR

**Normal Myocardium**

**Infarcted Myocardium**

**Ischemic Myocardium**

Contrast injection

First Pass

Delayed Enhancement

Rest --- Stress

time

SIEMENS63

# Myocardial Perfusion Imaging

pre-contrast

bolus injection

post-contrast

time

SIEMENS64

# First Pass Myocardial Perfusion

Saturation recovery TurboFLASH with parallel Imaging, GRAPPA rate 2, 1-4 slices/heart beat

SIEMENS65

# Ischemic Heart Disease – Ischemic Region versus Infarct

(a) Rest perfusion

(b) Stress perfusion

(c) Infarct

Sir Run Run Shaw Heart Center, Hong Kong

MR0862

SIEMENS66

# Current Pulse Sequences for Myocardial Perfusion with MR

- T1 weighted gradient echo sequences
- Steady State sequences
- Segmented echo planar sequences

SIEMENS
67

# Magnetization Preparation

- Inversion Recovery (180° pulse)
  - Very strong T1 contrast.
  - Multi-slice imaging difficult.
  - Sensitive to arrhythmia.
- Saturation Recovery (90° pulse)
  - Weaker T1 contrast.
  - Multi-slice capabilities.
  - Insensitive to arrhythmia.

SIEMENS
68

# Dynamic First-Pass Perfusion Saturation Recovery Preparation

- weaker T1 contrast
- Insensitivity to Arrhythmia
- Multi-slice capabilities
- CNR & SNR depend on TI

SIEMENS
69

# Outline

- Physiology and cardiac MR acquisitions
- Morphology
- Function
- Flow
- Perfusion
- Late enhancement

SIEMENS
70

# Late Enhancement Imaging

- Clinical goals
  - Identify locations of the myocardium that are 'dead' or injured, including infarcts, fibrous tissue due to myocarditis, intentional ablation, cardiomyopathies, etc
  - Identification from kinetic behavior of Gd-DTPA
- Sequence requirements
  - Highly T1 sensitive
  - High spatial resolution
  - Short acquisition time (obtain images in a BH)


SIEMENS
71

# Contrast Enhanced MR

SIEMENS
72

# Inversion Recovery Imaging


- Shorter T1 of contrast-enhanced infarcted region results in faster signal recovery following inversion.
- Nulling of signal from longer T1 viable myocardium results in maximum contrast ratio.



SIEMENS 75

## Effect of TI on Image Quality

- Imaging at the null-point of viable myocardium
- TI depends on:
  - Contrast dosage
  - Time between injection and imaging
  - Relaxation time between IR pulses (e.g., 2D versus 3D breathhold imaging)
- And TI is time dependent ...
  - Contrast in myocardium washes out with time






SIEMENS 77

## Inversion Time Scout

■ Inversion Time Scout Sequence

- Segmented IR-SSFP cine
- Quickly identify correct inversion time (TI).
- May be useful for T1 mapping





SIEMENS 81

### Summary

- Physiology and cardiac MR acquisitions
- Morphology
- Function
- Flow
- Perfusion
- Late enhancement

SIEMENS 82

### And for slide contributions....

### Thank you!

- Lon Simonetti
- Stefan Fischer
- Peter Kellman
- Yiu-Cho Chung
- Gerhard Laub
- Gary McNeal
- Kevin Johnson
- Carmel Hayes
- Dudley Pennell
- CH Luk