


### Measuring Brain Volume Changes: The Tools

Stephen Smith


Oxford University Centre for Functional MRI of the Brain  
(FMRIB)



Measuring Brain Volume Changes: The Tools : 1 of 31

### Measuring Brain Volume Changes: The Tools


- Tissue-type segmentation & bias field correction
- Longitudinal: temporal brain change
- Cross-sectional: single time point brain state
- Localised analyses




Measuring Brain Volume Changes: The Tools : 2 of 31

### Segmentation example: FAST from FSL

- First use BET to remove non-brain




- Input can be single image (e.g. T1, T2, Proton Density)
- Or several of these ("multi-channel")
- For multi-channel, all must be aligned (e.g. with FLIRT)




Measuring Brain Volume Changes: The Tools : 3 of 31

### MR Images - RF Inhomogeneity ("Bias Field")




- RF Inhomogeneity causes intensity variations in image
- Causes problems for simple threshold-based segmentation
- Need to remove bias field before or within segmentation



Measuring Brain Volume Changes: The Tools : 4 of 31

## Histograms - Tissue Intensity Distributions


- Histogram = "voxel count vs intensity"
- Mixture of Gaussians
- Model class means and widths
- If well separated, clear peaks, i.e. segmentation easy
- But overlap worsened by: bias, blurring, low resolution, head motion



Measuring Brain Volume Changes: The Tools : 5 of 31

## Initial Segmentation


- Need to bootstrap method somehow - initial segmentation
- Use "Tree-Structure K Means"
- Start with one class (Gaussian)
- Split & fit
- Re-split until enough classes
- Many other segmentations **only** do this step!



Measuring Brain Volume Changes: The Tools : 6 of 31

## Using Spatial Neighbourhood Information (MRF)


- Neighbourhood information: "if my neighbours are grey matter then I probably am too"
- Most methods (like the k-means initialisation) don't use spatial neighbourhood information
- Reduces noise and increases robustness
- Carried out using MRF (Markov Random Field) model



Measuring Brain Volume Changes: The Tools : 7 of 31

## At Last - the Overview!


- Initial segmentation - tree-structure k-means
- Iterate
  - Estimate bias field
  - Iterate segmentation
    - Update segmentation
    - Update tissue class parameters (mean and standard deviation)




Measuring Brain Volume Changes: The Tools : 8 of 31

## Examples - Single Channel Segmentation


Original,  
original plus bias field,  
histogram,  
best thresholding  
segmentation



No MRF used:  
estimated bias,  
segmentation, restored  
image, histogram




With MRF



Measuring Brain Volume Changes: The Tools : 9 of 31

## A-Priori Tissue Probability Maps

- A-priori maps created by averaging many aligned segmentations; can be used as priors in segmentation but can skew results
- If bias is very bad, priors can aid initial segmentation
- A-priori maps can also be optionally used to feed into final posteriors (e.g. to aid segmentation of deep gray); FAST doesn't use this by default, SPM2 does
- SPM5 improves the use of priors (compared with SPM2) by combining segmentation with alignment to priors.




Measuring Brain Volume Changes: The Tools : 10 of 31

## Longitudinal Change Analysis



Example: SIENA (Structural Image Evaluation, using  
Normalisation, of Atrophy, from FSL)

- Measures atrophy / general brain change
- Accurate and fully automatic
- Proven for a range of slice thicknesses
- Proven for a range of MRI sequences
- Accuracy 0.2% of brain volume
- Correction for scanner geometry drifts



Measuring Brain Volume Changes: The Tools : 11 of 31


## Example of Atrophy in Action



Measuring Brain Volume Changes: The Tools : 12 of 31

## BET : Brain Extraction Tool




- Eliminates all non-brain tissue
- Accurate and fully automatic
- < 1 minute processing time



Measuring Brain Volume Changes: The Tools : 13 of 31

## BET : Brain Extraction Tool


- Estimation of exterior surface of skull
- Used to hold image scale constant in registration



Measuring Brain Volume Changes: The Tools : 14 of 31

## FLIRT Linear Registration





- Register brains (full affine)
- Apply to skulls then re-register skulls, altering only scale and skew (corrects for changes in scanner geometry etc.)
- Apply to brains then re-optimize rotation and translation
- Use transform's midway position
- < 10 minutes for 3-step registration



Measuring Brain Volume Changes: The Tools : 15 of 31

## Atrophy Measurement Using Edge Motion


Find brain/non-brain edge points in image 1 using FAST tissue segmentation (including bias-field correction)



Measuring Brain Volume Changes: The Tools : 16 of 31

## Atrophy Measurement Using Edge Motion


At each edge point take 1D perpendicular profile



Measuring Brain Volume Changes: The Tools : 17 of 31

## Atrophy Measurement Using Edge Motion


- Find profile from same point in image 2
- Take windowed derivative of both
- Correlate with subvoxel accuracy



Measuring Brain Volume Changes: The Tools : 18 of 31

## Atrophy Measurement Using Edge Motion


- Thus for each edge point in image 1 the perpendicular motion is found
- Insensitive to changes in imaging
- Mean motion over whole brain surface calculated
- Convert into % brain volume change (PBVC) by estimating brain surface area and volume and then doing self-calibration
- Example: blue=atrophy, red="growth"



Measuring Brain Volume Changes: The Tools : 19 of 31

## Error Plots from "Normals"

- 16 normals, each scanned twice
- Range of slice thicknesses
- Error not dependent on slice thickness
- Error approximately 0.2%




Measuring Brain Volume Changes: The Tools : 20 of 31

## Error Plots from 3 Time Points from Patients

Data courtesy of V. Stevenson, D. Miller *et al.*, ION, London


- 39 patients (white-matter atrophy)
- Three time points
- Test accuracy by  
t0->t2  
vs  
t0->t1 + t1->t2
- Errors within 0.2%



Measuring Brain Volume Changes: The Tools : 21 of 31

## Cross-Sectional Atrophy: Brain State

- Example: SIENAX (from FSL)
- Measures brain volume normalised for head size
- Proven for T1, T2, PD etc.
- Test-retest error 0.5-1% (0.2-0.4% without normalisation)
- SE for normals 0.7% (1.5% without normalisation)




Measuring Brain Volume Changes: The Tools : 22 of 31

## SIENAX : Cross-Sectional Atrophy

- BET: find brain and skull
- FLIRT: use brain & skull to normalise to standard space
- Use standard space mask to cut brain stem and ensure no optic nerve/eyeballs
- Tissue segmentation (FAST, including partial volume estimation)
- Gives (normalised) brain (and grey & white) volume

Measuring Brain Volume Changes: The Tools : 23 of 31

## Normalised Brain Volume vs DoB (Normals)



Measuring Brain Volume Changes: The Tools : 24 of 31



## SIENAX Results : MS vs Controls

Group	N	Mean Volume /10 <sup>6</sup> mm <sup>3</sup>	SD /10 <sup>6</sup> mm <sup>3</sup>	P
Controls	20	1.45	0.05	
Total MS	72	1.39	0.10	0.0001
EDSS <2	33	1.44	0.07	
EDSS <5	61	1.41	0.09	0.01
EDSS 2-4	27	1.38	0.10	0.0001
EDSS 5-8	10	1.27	0.10	0.0001

- Data from Siena, Italy
- ANOVA Tukey shows significant difference between controls and all EDSS bins except for EDSS<2

Measuring Brain Volume Changes: The Tools : 25 of 31

## Voxelwise Cross-Subject Statistics - SIENAr

Extend SIENA for voxelwise cross-subject statistics, e.g.:

- where is atrophy different between two groups,
- or where does atrophy correlate with EDSS
- run SIENA to get edge "flow" image
- dilate
- transform to standard space
- mask with standard-space edge-mask
- blur (optional)
- carry out voxelwise cross-subject statistics




Measuring Brain Volume Changes: The Tools : 27 of 31

## Cross-Subject Comparisons (Partial Head Images)


If different subjects have varying amounts of head in field of view, atrophy comparisons are not valid. Therefore use standard-space-based masking or Z limits so that all is consistent.

For example, to only use data where  $-20\text{mm} < Z < +50\text{mm}$  (in standard space), use options `-b -20 -t 50`



Measuring Brain Volume Changes: The Tools : 26 of 31

## Voxelwise Cross-Subject Statistics - SIENAr




Example - one-group MS atrophy

Measuring Brain Volume Changes: The Tools : 28 of 31

## Regional Atrophy Measurement - SIENAX

With SIENAX, carry out regional breakdown

This gives peripheral GM volume and ventricular CSF volume, using standard-space masks.



Measuring Brain Volume Changes: The Tools : 29 of 31

## Voxelwise Volumetry - VBM

- VBM: Voxelwise structural stats, e.g. with SPM or SIENAX+IRTK(nonlinear reg):
- Segment -> GM PVE
- Align to standard space (with density modulation)
- Voxelwise cross-subject stats
- Pros: fully automated, easy to test whole brain
- Con: ambiguity between cross-subject geometry shifts and intensity changes...hard to interpret results sometimes

Measuring Brain Volume Changes: The Tools : 30 of 31

## References

- Y. Zhang, M. Brady, and S. Smith.  
Segmentation of brain MR images through a hidden Markov random field model and the expectation maximization algorithm.  
*IEEE Trans. on Medical Imaging*, 20(1):45-57, 2001.
- S.M. Smith.  
Fast robust automated brain extraction.  
*Human Brain Mapping*, 17(3):143-155, November 2002.
- S.M. Smith, Y. Zhang, M. Jenkinson, J. Chen, P.M. Matthews, A. Federico, and N. De Stefano.  
Accurate, robust and automated longitudinal and cross-sectional brain change analysis.  
*NeuroImage*, 17(1):479-489, 2002.
- M. Jenkinson, P.R. Bannister, J.M. Brady, and S.M. Smith.  
Improved optimisation for the robust and accurate linear registration and motion correction of brain images.  
*NeuroImage*, 17(2):825-841, 2002.
- S.M. Smith, M. Jenkinson, M.W. Woolrich, C.F. Beckmann, T.E.J. Behrens, H. Johansen-Berg, P.R. Bannister, M. De Luca, I. Drobnjak, D.E. Flitney, R. Niazy, J. Saunders, J. Vickers, Y. Zhang, N. De Stefano, J.M. Brady, and P.M. Matthews.  
Advances in functional and structural MR image analysis and implementation as FSL.  
*NeuroImage*, 23(S1):208-219, 2004.

Measuring Brain Volume Changes: The Tools : 31 of 31